



# CICV\* ? Use Ventrain !

full ventilation by active expiration:

- » MV >6 l/min
- » CO<sub>2</sub> removal

through a 2 mm lumen  
even with an obstructed upper airway



In case of a CICV (cannot intubate cannot ventilate) situation Ventrain can be connected to a 2mm transtracheal catheter.

**Ventrain starts a new era in ventilation**

[www.ventrain.com](http://www.ventrain.com)

Another quality product  
selected and stocked by:  
W.M. Bamford & Co Ltd  
p 0800 5739109  
e [wayne@bamford.co.nz](mailto:wayne@bamford.co.nz)  
[www.bamford.co.nz](http://www.bamford.co.nz)



when performance counts Bamford delivers