

Established 1974

Austmel

assuring your future

WASH-CHECKS™ PRO

Rapid test to detect protein on surfaces

Re-order Code: WC112


- **WASH-CHECKS PRO HAS THE ABILITY TO DETECT RESIDUAL PROTEINS LEFT BEHIND ON THE SURFACES OF WASHER-DISINFECTORS, ULTRASONIC CLEANERS, ENDOSCOPES AND SURGICAL INSTRUMENTS**

Can detect protein residues within 1 µg sensitivity.

If protein is present, the reagent will turn blue within 10 seconds.

The higher level of contamination of protein, the darker and faster the colour changes to blue.

- **SIMPLE, FAST, EASY TO USE**
- **COST EFFECTIVE**



ORDERING INFORMATION

PART #	PACKAGING
WC112	25 tests/kit 2 kits/case

T 1800 074 085 Customer Support

