

List of Disinfectants



For disinfecting our products from the product lines



we recommend the following disinfecting agents:

Name	Manufacturer	Active Ingredient	Listed in List of VAH*	Concentration and Contact Time
Cidex OPA	Ethicon GmbH	aldehyde	yes	For using the individual disinfecting agents, please pay attention to the respective manufacturer's data!
Gigasept FF	Schülke & Mayer GmbH	aldehyde, separation of aldehyde	yes	
Korsolex basic	Bode Chemie GmbH & Co.	aldehyde separation of aldehyde	yes	
Korsolex extra	Bode Chemie GmbH & Co.	quaternary compounds, aldehyde, separation of aldehyde	yes	
Korsolex PAA	Bode Chemie GmbH & Co.	hydrogen peroxide and peracetic acid	no	
Sekusept aktiv	Ecolab Deutschland GmbH	peroxide compounds	yes	
Sekusept pulver classic	Ecolab Deutschland GmbH	peroxide compounds	yes	
Sekusept plus	Ecolab Deutschland GmbH	alkylamine derivative	yes	
Sekusept forte	Ecolab Deutschland GmbH	quaternary compounds, aldehyde	yes	

After each disinfectant process and before each use our medical products must be checked regarding integrity and functional efficiency.

* List of disinfectants of the "Verbund für angewandte Hygiene e.V." (Association for Applied Hygiene)