



Lohmann & Rauscher

Compression Catalogue

Complete range of compression solutions



Rosidal®
Short stretch
bandage range



ReadyWrap™
Adjustable
compression wraps



VENOSAN®
Compression
stockings



Table of Contents

Compression Bandaging

Rosidal® TCS	5
Rosidal® sys	6
Rosidal® K	7
Rosidal® SC	8
Rosidal® CC	9
Rosidal® soft	10
Rosidal® haft	11
Cellona® Synthetic Undercast Padding	12
tg® grip	13
tg® soft	14
Sama™ Frottee Tubular Padding	15
Komprex™	16

Compression Stockings and Garments

Product lines	19
Product overview	20
VENOSAN® 3000 COTTON	22
VENOSAN® 4000	23
VENOSAN® 5000	24
VENOSAN® 6000	25
VENOSAN® 7000	26
VENOSAN® 8000 ULCERFIT	27
VENOSAN® 4000 Armsleeves	28
VENOSAN® 7000 Armsleeves	29
VENOSAN® Flat	30

Support Stockings

Product overview	34
VENOSAN® LEGLINE 20	36
VENOSAN® LEGLINE 30	37
VENOSAN® MICROFIBERLINE	38
VENOSAN® SILVERLINE	39
VENOSAN® ATHLETIC	40
VENOSAN® SPORTSLINE	41

Anti-Embolic Stockings

VENOSAN® AES	43
--------------	----

Adjustable Compression Wraps

ReadyWrap™	45
TributeWrap™	46

Additional Compression Garments

Caresia™	48
ExoSoft™	49
ExoStrong™	50
SwellSpots™	51

Accessories for compression stockings

VENO-GLIDER	53
VENOSAN® Dot Gloves	53
COLTEX	53
VENOSAN® Eureka On	53

Measuring and ordering information

VENOSAN® Measuring Guide	55
VENOSAN® Sizing Guide	56
VENOSAN® Measuring and Sizing Guide for Arm Compression Stockings	57
VENOSAN® 3000 COTTON	58
VENOSAN® 4000	59
VENOSAN® 5000	62
VENOSAN® 6000	64
VENOSAN® 7000	66
VENOSAN® 8000 ULCERFIT	68
VENOSAN® 4000 Armsleeves	69
VENOSAN® 7000 Armsleeves	70
VENOSAN® LEGLINE 20	71
VENOSAN® LEGLINE 30	72
VENOSAN® MICROFIBERLINE	73
VENOSAN® SILVERLINE	74
VENOSAN® ATHLETIC	75
VENOSAN® SPORTSLINE	76
VENOSAN® AES	77
ReadyWrap™, TributeWrap™, Caresia™, Exo, Measuring Instructions	78
ReadyWrap™ Upper Extremity Garments	79
ReadyWrap™ Lower Extremity Garments	80
TributeWrap™ Upper & Lower Extremity Garments	83
Caresia™ Upper Extremity Garments	84
Caresia™ Lower Extremity Garments	85
SwellSpots™	86
ExoSoft™	87
ExoStrong™ Upper Extremity Garments	88
ExoStrong™ Lower Extremity Garments	89



Lohmann & Rauscher

Compression Bandaging





Return to contents page

Rosidal® TCS

Two component compression system



Product benefits

- Rosidal TCS is made up of:
 - soft compression bandage (Rosidal SC) = first layer
 - cohesive compression bandage (Rosidal CC) = second layer
- Two components combined in a handy, ready-for-use set
- Integrated skin protection in the reusable first layer
- Short-stretch cohesive compression layer provides effective therapeutic compression
- Cohesive quality of the bandage allows for extended rebandaging intervals without slippage for up to seven days
- Offers greater comfort at rest which improves patient concordance and encourages mobility and confidence
- Safe-Loc: safe to apply at full stretch and no pressure indicators required, ensuring consistent and accurate application
- Simple and easy to apply
- Thin and supple, allowing good mobility in the ankle joint and shoes can be worn comfortably
- Latex-free
- Extensive clinical evidence

Indications

- For the treatment of lower leg ulcers and concomitant conditions
- For the treatment of chronic oedema



Rosidal® TCS

System components	Size	REF	Shipping units (pcs.)
Component 1: Rosidal SC	10 cm x 3.5 m	26484	1 / 14
Component 2: Rosidal CC	10 cm x 6.0 m		



Return to contents page

Rosidal® sys

Short stretch compression system



Product benefits

- Effective and economical as padding and compression bandages can be reused. (Rosidal soft and Rosidal K)
- Ready-for-use products combined in one handy set
- The system can be worn continuously for up to 7 days (high level of compression pressure is maintained throughout)
- Components: tg Tubular, Rosidal Soft (padding), Rosidal K (short stretch compression), Porofix (adhesive tape), Mollelast haft (cohesive conforming bandage)

Indications

- For the treatment of lower leg ulcers



Rosidal® sys

Number	Contents	Contents Size	REF	Shipping units (pcs.)
1	tg® Tubular Bandage	Size 7	16984	1 / 4
4	Rosidal® soft	10 cm x 2.5 m x 0.2 cm		
2	Rosidal® K	8 cm x 5 m		
2	Rosidal® K	10 cm x 5 m		
1	Porofix®	2.5 cm x 5 m		
1	Mollelast® haft	8 cm x 20 m		



Return to
contents page

Rosidal® K

Short stretch bandage



Product benefits

- 100% cotton short stretch bandage
- Reusable; washable up to 50 times
- Safe-Loc: safe to apply at full stretch and no pressure indicators required
- Non-cohesive
- Easy application
- Air-permeable

Indications

- For powerful compression (e.g. in phlebology)
- To immobilise the limbs
- For support and pressure relief in traumatology and sports medicine
- To retain dressings



Rosidal® K

Size	REF	Shipping units (pcs.)
5 m (stretched length), with bandage clips (except for 4 cm), individually packed in a folding box		
4 cm x 5 m	22199	1 / 125
6 cm x 5 m	22200	1 / 108
8 cm x 5 m	22201	1 / 108
10 cm x 5 m	22202	1 / 90
12 cm x 5 m	22203	1 / 54
10 m (stretched length), without bandage clips, individually packed in a folding box		
12 cm x 10 m	22205	1 / 40
10 cm x 10 m	22206	1 / 60



Return to contents page

Rosidal® SC

Soft compression bandage



Product benefits

- Helps distribute pressure evenly
- Integrated skin protection layer which provides:
 - Improved comfort and breathability
 - Reduced skin irritation
 - No need for an additional tubular bandage
- Light compression effect
- Reusable: washable up to 50 times
- First layer of Rosidal TCS
- Permeable to air and water vapour

Indications

- Can be used as the first compression layer to protect the skin and to make the bandage less prone to slipping
- For light compression under lymphological or phlebological compression bandages



Rosidal® SC

with bandage clips, individually packed in a folding box

Size	REF	Shipping units (pcs.)
10 cm x 2.5 m	33500	1 / 18
15 cm x 2.5 m	33501	1 / 24
10 cm x 3.5 m	33503	1 / 24



Return to contents page

Rosidal® CC

Cohesive compression bandage



Product benefits

- High static stiffness
- Delivers high working, low resting pressure
- Latex-free
- Cohesive on both sides
- Safe-Loc: safe to apply at full stretch and no pressure indicators required
- Simple and easy to apply
- Top layer of Rosidal TCS

Indications

- To immobilise the limbs
- For the retention of multi-layer bandages
- For compression



Rosidal® CC

6 m (stretched length), individually packed in a cardboard box

Size	REF	Shipping units (pcs.)
10 cm x 6 m	33499	1 / 24



Return to contents page

Rosidal® soft Foam bandage



Product benefits

- Lengthwise and widthwise elastic
- Ensures even distribution of pressure
- Increases patient comfort (e.g. cushioning of the tibial crest)
- Reduces slippage
- Highly permeable to air and humidity
- Tear-resistant and stretchable
- Reusable: washable up to 50 times
- Maintains its elasticity and cushioning effect even after washing

Indications

- For padding under lymphological and phlebological compression bandages



Rosidal® soft

2.5 m (unstretched length), with bandage clips, individually packed in folding boxes

Thickness	Dimensions	REF	Shipping units (pcs.)
0.3 cm	10 cm x 2.5 m x 0.3 cm	23110	1 / 30
0.4 cm	10 cm x 2.5 m x 0.4 cm	23111	1 / 24
0.4 cm	12 cm x 2.5 m x 0.4 cm	23112	1 / 16
0.4 cm	15 cm x 2.5 m x 0.4 cm	23113	1 / 14

Rosidal® soft

2.0 m unstretched length, with bandage clips, packed in pairs in a folding box

Strength (cm)	Dimensions (cm x m x cm)	REF	Shipping units (pcs.)
0.2 cm	10 cm x 2 m x 0.2 cm	13016	1 / 23



Return to contents page

Rosidal® haft

Cohesive short stretch compression bandage



Product benefits

- Delivers high working, low resting pressure
- Stays reliably in place
- Compression force is well tolerated due to the high degree of static stiffness
- High compression is sustained over a longer period
- Cohesive on both sides
- Safe-Loc: safe to apply at full stretch and no pressure indicators required
- Simple and easy to apply

Indications

- For compression of the limbs for phlebological and lymphatic conditions
- For the treatment of oedema
- After sclerotherapy and varicose vein surgery
- For support and pressure relief of the capsular ligament system



Rosidal® haft

5 m (stretched length), individually packed in a cardboard box non-sterile

Width	REF	Shipping units (pcs.)
6 cm x 5 m	31973	1 / 48
8 cm x 5 m	31974	1 / 32
10 cm x 5 m	31975	1 / 42
12 cm x 5 m	31976	1 / 24



Return to contents page

Cellona® Synthetic undercast padding



Versatile and protective

Cellona Synthetic Undercast Padding can be used for a number of different applications. In addition to its protective function under plaster of Paris and synthetic casts, it is used as padding for splints and compression bandages.

Flexibility and comfort

The fine, fibrous web enables smooth, spiral application but can also be teased apart where required, preventing bulging, to enhance wearing comfort.

The material (100% polyester fibres) is well suited to be used in narrow areas as it can be tapered by gently stretching it lengthwise.

An adhesive-like effect inherent to the material makes the padding easier to apply and smooth out.



Cellona Synthetic Undercast Padding

Size	REF	Shipping units (pcs.)
5 cm x 2.7 m	34591	6 / 144
7.5 cm x 2.7 m	34592	6 / 72
10 cm x 2.7 m	34593	6 / 48
15 cm x 2.7 m	34594	6 / 72



Return to contents page

tg[®] grip Tubular support bandage



Product benefits


- Supple material with uniform compression effect
- Air-permeable
- Light padding effect
- Application with virtually no creases
- Knitted seamlessly
- Widthwise stretch
- Skin-friendly due to high percentage of cotton
- Washable up to 60°C and sterilisable by steam (134°C)
- Compression can be enhanced with a second layer of material or by using a smaller size

Product composition

- 85% cotton, elastodiene (latex); covered with a polyamide yarn

Indications

- For support and retention with a light compression effect



tg[®] grip

10 m roll in a cardboard dispenser

	Size	Description	REF	Shipping units (pcs.)
A	4.6 cm	children's hands, children's forearms	24120	1 / 24
B	6 cm	children's arms, children's lower legs	24121	1 / 16
C	6.5 cm	slim arms, children's legs	24122	1 / 16
D	7.5 cm	medium-sized arms, lower legs	24123	1 / 16
E	8.75 cm	large arms, legs, slim thighs	24124	1 / 12
F	10 cm	legs, medium-sized thighs	24125	1 / 12
G	12 cm	legs, for large thighs	24126	1 / 8
J	17.5 cm	small trunks	24127	1 / 6
K	21.5 cm	medium-sized trunks	24128	1 / 4
L	32.5 cm	large trunks	24129	1 / 4



Return to contents page

tg[®] soft Tubular padding bandage



Product benefits


- Soft and very elastic
- Easy to cut to the required size
- Virtually no creasing when applied, skin-friendly
- Conforms to any body contour
- Adheres to the cast
- Highly permeable to water vapour and air
- Low absorption capacity ensures high wearing comfort
- Padding and skin protection in one

Product composition

- 85% cotton (unbleached), 15% elastane

Indications

- For padding under plaster and synthetic casts



tg[®] soft

10 m roll in a cardboard dispenser

Size	Use	REF	Shipping units (pcs.)
S	forearm, child's arm, child's leg	13200	1 / 8
M	lower leg, arm	13201	1 / 6
L	leg	13202	1 / 4



Return to contents page

Sama™ Frottee Tubular padding



Product benefits

- Easy to apply
- Soft
- Increased comfort from padding
- Will not crease or wrinkle
- Allows skin to breathe
- Reusable - washable up to 20 times
- Contours to all shapes, sizes and lengths

Product composition

- 95% cotton, 3% polyamide, 2% elastane

Sama™ Frottee Tubular Padding

10 m roll in a cardboard dispenser

Colour	Size	Circumference	REF	Shipping units (pcs.)
White	Small	15 - 35 cm	159366	1
White	Medium	30 - 55 cm	159367	1
White	Large	40 - 75 cm	159368	1
White	X-Large	50 - 95 cm	159369	1
Pink	Small	15 - 35 cm	159370	1
Pink	Medium	30 - 55 cm	159371	1
Pink	Large	40 - 75 cm	159372	1
Blue	Small	15 - 35 cm	159374	1
Blue	Medium	30 - 55 cm	159375	1
Blue	Large	40 - 75 cm	159376	1
Green	Small	15 - 35 cm	159378	1
Green	Medium	30 - 55 cm	159379	1
Green	Large	40 - 75 cm	159380	1



Return to
contents page

Komprex™ Foam pads & bandages



Product benefits

Komprex

- Can be cut to size individually
- Permeable to air and water vapour

Komprex II

- Can be cut to size individually
- Breathable and non-slip
- Special wavy surface structure
- Micro-massaging effect particularly benefits patients with lymphoedema

Indications

Komprex

- Komprex Foam Rubber Bandage
– for compression in oedema treatment in phlebology and lymphology
- Komprex Foam Rubber Pads
– for support with functional bandages in sports medicine
- Use in combination with other compression products to change local pressure

Komprex II

- Padding with compression bandages

Komprex™ Foam Rubber Pad

individually packed in a folding box

Size	Thickness	Dimensions	REF	Shipping units (pcs.)
0 kidney-shaped	10 mm	9 cm x 5 cm	22301	1 / 75
1 kidney-shaped	12 mm	12 cm x 6.5 cm	22302	1 / 30

Komprex™ Foam Rubber Bandage

rolled, individually packed in a folding box

Thickness (mm)	Dimensions	REF	Shipping units (pcs.)
5 mm	8 cm x 1 cm x 0.5 cm	22310	1 / 30
10 mm	10 cm x 1 cm x 1 cm	22312	1 / 20
5 mm	8 cm x 2 cm x 0.5 cm	22313	1 / 30
10 mm	8 cm x 2 cm x 1 cm	22314	1 / 15

Komprex™ II Foam Rubber Sheet

wavy, individually packed in a foil pouch

Dimensions	REF	Shipping units (pcs.)
65 cm x 65 cm	16985	1 / 6



Lohmann & Rauscher

Compression Stockings and Garments



Simple. Effective. Swiss quality.



Return to
contents page

VENOSAN® Compression Stockings help to make life easier in many situations. The pressure VENOSAN Compression Stockings exert promotes healthy veins. They can help to prevent venous disease as well as being an important part of the treatment for patients being treated for acute and chronic venous disorders.

The diverse VENOSAN range offers the right product for a wide variety of needs. All VENOSAN Compression Stockings have the following characteristics in common:

Simple.

The VENOSAN range provides a simple guide to help you determine the correct size, compression class and style of compression stocking.

The wear and care instructions included with VENOSAN Compression Stockings are easy to follow which patients appreciate.

Effective.

VENOSAN Compression Stockings are RAL-certified, meaning that they provide effective therapeutic compression. The effective design of VENOSAN Compression Stockings means that for patients, many will experience enhanced treatment success and improved quality of life.

Swiss quality.

Products from Switzerland are known for excellent workmanship and top quality.

Manufactured using high quality materials, such as the original Lycra® and elastane, VENOSAN Compression Stockings offer wearing comfort and durability.



Product lines

Medical compression stockings and pantyhose.



Return to contents page

Line	Colours	CCL									
			AD Below knee stocking	ADH Below knee stocking with dotted top band	AG Thigh- length stocking	AGH Thigh-length stocking with dotted top band	AGH Thigh-length stocking with lace top band	AGG Thigh- length stocking with belt	AT Pantyhose	ATU Maternity pantyhose	
VENOSAN® 3000 NATURAL LOOK	Beige	CCL 2	●	●		●					
		CCL 3	●								
VENOSAN® 4000 CLIMATE EFFECT	Mexico	CCL 1	●	●	○	●	●	○	●	●	
	Marokko	CCL 2	●	●	○	●	●	●	●	●	
	Black	CCL 2	●	●	○	●	●	●	●	●	
VENOSAN® 5000 THE MAGIC OF SILVER	Mexico	CCL 1	●	●	○	●	●	●	●	○	
	Marokko	CCL 2	●	●	○	●	●	●	●	○	
	Black	CCL 3	●	○	○	○	○	○	○	○	
VENOSAN® 6000 CLASSIC	Beige	CCL 1	●			●			●		
	Black	CCL 2	●			●		●	●		
		CCL 3	●			●			●		
VENOSAN® 7000 NATURAL TOUCH	Beige	CCL 1	●								
	Black	CCL 2	●	●	○	●	○	●	○		
		CCL 3	●	●	○	●	○	○	○		
VENOSAN® 8000 ULCERFIT SUCCESSFUL THERAPY	Beige	CCL 2	●								
		CCL 3	●								
VENOSAN® SOFT VENOSAN® PRO VENOSAN® STRONG FLAT KNIT THERAPY	Beige	CCL 1	○	○	○	○		○	○		
	Black	CCL 2	○	○	○	○		○	○		
		CCL 3	○	○	○	○		○	○		
		CCL 4	○	○	○	○		○	○		

● available in standard sizes ○ available in custom sizes only ● available in standard and custom sizes

Pressure at ankle **CCL 1** 18–21mmHg **CCL 2** 23–32mmHg **CCL 3** 34–46mmHg **CCL 4** >49mmHg

Product overview

Medical compression stockings and pantyhose.



Return to contents page

VENOSAN® 3000 COTTON NATURAL LOOK

- Durable
- 100% cotton on the skin
- Absorbs moisture
- Easy to don

Indications

- CVI (chronic venous insufficiency) 2–3
- Lipoedema stage 2–3
- After vein surgery
- PTS (post-thrombotic syndrome)
- Mild lymphoedema



Composition


- CCL 2 49% cotton
27% nylon polyamide
24% elastane
- CCL 3 40% cotton
27% nylon polyamide
33% elastane

VENOSAN® 4000 CLIMATE EFFECT

- Climate effect
- Semi-transparent
- Comfortable
- Highly elastic

Indications

- CVI (chronic venous insufficiency) 1
- Varicosis during pregnancy
- After vein surgery



Composition

- 71% TACTEL® polyamide
- 29% elastane

VENOSAN® 5000 THE MAGIC OF SILVER

- Odour-inhibiting
- Anti-static
- Temperature-regulation
- Anti-bacterial

Indications

- CVI (chronic venous insufficiency) 1–2
- After vein surgery
- Mild skin conditions
- Lipoedema stage 1



Composition

- 40% TACTEL® polyamide
- 30% elastane
- 20% cotton
- 10% polyamide/silver (X-Static®)

Product overview

Medical compression stockings and pantyhose.



Return to contents page

VENOSAN® 6000 CLASSIC

- Maximum sheerness
- Breathable
- Ease of care

Indications

- CVI (chronic venous insufficiency) 1–2
- Varicose veins
- Varicosis during pregnancy
- After vein surgery



Composition


70% nylon polyamide
30% elastane

VENOSAN® 7000 NATURAL TOUCH

- High stiffness (resistance)
- Skin-friendly
- Temperature regulation

Indications

- CVI (chronic venous insufficiency) 2–3
- Lipoedema stage 2–3
- After vein surgery
- After thrombosis
- PTS (post-thrombotic syndrome)
- Mild lymphoedema
- Marked orthostatic stress



Composition

62% TACTEL® polyamide
26% elastane
8% cotton
4% Lyocell/SeaCell®

VENOSAN® 8000 ULCERFIT SUCCESSFUL THERAPY

- Two-layer stocking system
- Fast and easy to use
- High wearing comfort
- Liner can be worn at night

Indications

- Treatment of venous leg ulcers
- Helps to prevent recurrence of venous leg ulcers



Composition

Overstocking:
62% TACTEL® Polyamide
26% Elastane
8% Cotton
4% Lyocell/SeaCell®

Liner:
90% Polyamide
10% Elastane

VENOSAN® SOFT VENOSAN® PRO VENOSAN® STRONG FLAT KNIT THERAPY

- Flat-knit custom made garments
- Smooth and durable
- Low to high stiffness
- Manufactured to European standards

Indications

- Treatment and management of Chronic Oedema and Lymphoedema



Composition

Polyamide & Elastane

VENOSAN PRO (CCL II)
& VENOSAN STRONG (CCL 4S)
also contains cotton

VENOSAN PRO also contains Latex



Return to contents page

VENOSAN® 3000 COTTON

NATURAL LOOK

100% cotton - natural and soft.



■ Natural

Cotton allows a higher degree of wearer comfort for all those who prefer a natural fibre against the skin, or who are allergic to synthetic fibres. 100% cotton against the skin.

■ Durable

The construction ensures the compression stocking is durable and hard wearing making it ideal for long-term compression therapy.

■ Easy to don

The stretchable knitted cotton fabric allows easy donning of the compression stocking.

[Click for VENOSAN 3000 ordering and sizing details](#)

Compression Stockings and Garments

Beige



CCL 2

CCL 3

Garment Styles



Indications

- CVI 2–3 (chronic venous insufficiency)
- Lipoedema stage 2–3
- Post vein surgery
- PTS (post-thrombotic syndrome)
- Mild lymphoedema

Composition

- CCL 2: 49% cotton, 27% nylon polyamide, 24% elastane
- CCL 3: 40% cotton, 27% nylon polyamide, 33% elastane



Return to contents page

VENOSAN® 4000

CLIMATE EFFECT

Microfibre - fine and stylish.



■ Climate effect

The air permeable yarns and knitting construction of this stocking allow for moisture to be wicked away from the skin and through the stocking.

■ Fashionably elegant

Wearable in all seasons. It has a contemporary stylish look, with the sheer transparency.



■ Comfortable

The velvety touch of the high tech TACTEL® microfibrils offers outstanding comfort. The microfibrils are many times thinner than a human hair.

■ Made to Measure available

[Click for VENOSAN 4000 ordering and sizing details](#)

Compression Stockings and Garments



Garment Styles



Indications

- CVI 1 (chronic venous insufficiency)
- Pregnancy-induced varicose veins
- Milder symptoms
- Post vein surgery

Composition

- 71% TACTEL® polyamide
- 29% elastane



Return to contents page

VENOSAN® 5000

THE MAGIC OF SILVER

For ideal temperature regulation.



■ Durable

The innovative yarn combination of silver and cotton make this stocking suitable for active men and women.

■ Regulates temperature

Silver is well known for its excellent heat conductive properties. In combination with the microfibre TACTEL®, silver produces rapid dissipation of skin moisture, helping to ensure optimal wearing comfort.



■ Silver-coated X-Static® fibres

Contains silver-coated X-Static fibre. No reduction of the properties of silver even after 150 washing cycles.

■ Made to Measure available

[Click for VENOSAN 5000 ordering and sizing details](#)

Compression Stockings and Garments



Garment Styles



Indications

- CVI 1–2 (chronic venous insufficiency)
- Post vein surgery
- Mild skin conditions
- Lipoedema stage 1

Composition

- 40% TACTEL® polyamide
- 30% elastane
- 20% cotton
- 10% polyamide / silver X-Static®



Return to contents page

VENOSAN® 6000

CLASSIC



Reliable and cost-effective.



- **Durable knit construction**
Basic yarns make the stocking cost-effective.
- **Anatomically knitted**
For maximum comfort, does not slip.
- **Ease of care**
Can be machine washed at 40°C on a delicates programme. Dries naturally overnight.

[Click for VENOSAN 6000 ordering and sizing details](#)

Compression Stockings and Garments



Garment Styles



Indications

- CVI 1-2 (chronic venous insufficiency)
- Varicose veins
- Varicosis during pregnancy
- After vein surgery

Composition

- 70% nylon polyamide
- 30% elastane



Return to contents page

VENOSAN® 7000

NATURAL TOUCH

High stiffness with optimal comfort.



■ High stiffness

Due to the enhanced stiffness the stockings provide a massage effect and are characterised as shaping, opaque compression stockings with high stiffness.

■ Skin-friendly

The composition ensures both high resistance and optimal wearing comfort.


■ Strong natural fibres

Due to the durability of the material these stockings are suitable for people who demand a lot from their stockings at work, in sports and in their leisure time.

■ Made to Measure available

[Click for VENOSAN 7000 ordering and sizing details](#)

Compression Stockings and Garments



Garment Styles



Indications

- CVI 2–3 (chronic venous insufficiency)
- Lipoedema stage 2–3
- Following acute thrombosis (following dressing)
- Post vein surgery
- PTS (post-thrombotic syndrome)
- Mild lymphoedema
- Marked orthostatic stress

Composition

- 62% TACTEL® polyamide
- 26% elastane
- 8% cotton
- 4% Lyocell/SeaCell® (seaweed fibre)



Return to contents page

VENOSAN® 8000 ULCERFIT

SUCCESSFUL THERAPY

Compression Stocking System.



- **Two layer system**
Offers optimum therapy thanks to the ideal combination of liner and overstocking.
- **Time saving treatment**
Fast and easy to use, patient can apply the system without assistance.
- **Better patient compliance**
High wearer comfort. The flat and thin surface of the liner ensures easy donning and doffing of the overstocking.

[Click for VENOSAN 8000 ordering and sizing details](#)

Compression Stockings and Garments

Beige



Version



Indications

- Treatment of venous leg ulcers
- Helps to prevent recurrence of venous leg ulcers

Composition

Overstocking:

- 62% TACTEL® Polyamide
- 26% Elastane
- 8% Cotton
- 4% Lyocell/SeaCell®

Liner:

- 90% Polyamide
- 10% Elastane





Return to contents page

VENOSAN® 4000 Armsleeves

CLIMATE EFFECT

Microfibre - fine and stylish.



■ Climate effect

The air permeable yarns and knitting construction of this stocking allow for moisture to be wicked away from the skin and through the stocking.

■ Fashionably elegant

Wearable in all seasons. It has a contemporary stylish look, with the sheer transparency.

■ Comfortable


The velvety touch of the high tech TACTEL® microfibrils offers outstanding comfort. The microfibrils are many times thinner than a human hair.

■ Made to Measure available

[Click for VENOSAN 4000 Armsleeves ordering and sizing details](#)

Compression Stockings and Garments

Beige



Garment Styles



Indications

- Mild lymphoedema

Composition

- 71% TACTEL® polyamide
- 29% elastane



Return to contents page

VENOSAN® 7000 Armsleeves

NATURAL TOUCH

High stiffness with optimal comfort.



■ High stiffness

Due to the enhanced stiffness the stockings provide a massage effect and are characterised as shaping, opaque compression stockings with high stiffness.

■ Skin-friendly

The composition ensures both high resistance and optimal wearing comfort.

■ Strong natural fibres

Due to the durability of the material these stockings are suitable for people who demand a lot from their stockings at work, in sports and in their leisure time.

■ Made to Measure available


[Click for VENOSAN 7000 Armsleeves ordering and sizing details](#)

Compression Stockings and Garments

Beige



Garment Styles



Indications

- Mild lymphoedema

Composition

- 62% TACTEL® polyamide
- 26% elastane
- 8% cotton
- 4 % Lyocell/SeaCell® (seaweed fibre)



Return to contents page

VENOSAN® Flat

Flat knit compression garments Made To Measure



Made-to-measure flat knit compression garments for the treatment of Lymphoedema and Chronic Oedema.

- Extremely high quality, made-to-measure flat knit compression garments.
- Low to high stiffness.
- Manufactured to European standards.
- RAL certified (indent).
- Oeko-Tex Standard 100 (indent).
- Comprehensive range of upper and lower extremity stockings.
- 3 types of fabrics – enhanced ability to select the fabric that best suits your patients' needs.

VENOSAN Soft

VENOSAN Pro

VENOSAN Strong



Most common garment styles.

See following page for complete range



Indications

- Treatment and management of chronic oedema and lymphoedema

Composition

Polyamide & Elastane
VENOSAN PRO (CCL II)
& VENOSAN STRONG (CCL 4S)
also contains cotton
VENOSAN PRO also contains Latex

VENOSAN® Flat

Range overview



Return to
contents page

		VENOSAN® SOFT <small>Latex Free</small>			VENOSAN® PRO		VENOSAN® STRONG <small>Latex Free</small>		
Colour range		● Beige	● Black		● Beige	● Black		● Beige	
Compression Class		Class 1 (18-21 mmHg)	Class 2 (23-32 mmHg)	Class 3 (34-46 mmHg)	Class 2 (23-32 mmHg)	Class 3 (34-46 mmHg)	Class 3 (34-46 mmHg)	Class 4 (>49 mmHg)	Class 4S (>49 mmHg)
Lower Extremity									
AB	Ankle	●	●	●	●	●	●	-	-
AD	Below Knee	●	●	●	●	●	●	●	●
AF	Mid-Thigh	●	●	●	●	●	●	●	●
AG	Thigh High	●	●	●	●	●	●	●	●
AGT	Thigh High w. Waist Attachment	●	●	●	●	●	●	●	●
AT	Pantyhose	●	●	●	●	●	●	●	●
ATS	Single Leg Pantyhose, short side finishing at point F or G	●	●	●	●	●	●	●	●
BD	Ankle To Knee	●	●	●	●	●	●	●	●
BT	Ankle To Waist	●	●	●	●	●	●	●	●
DF	Knee	●	●	●	●	●	●	●	●
DG	Knee To Thigh	●	●	●	●	●	●	●	●
DGT	Knee To Thigh w. Hip Attachment	●	●	●	●	●	●	●	●
DT	Knee To Waist (Capri)	●	●	●	●	●	●	●	●
DTS	Single Leg Capri, short side finishing at point F or G	●	●	●	●	●	●	●	●
FT	Bermuda Pantyhose	●	●	●	●	●	●	●	●
FTGT	Single Leg Bermuda Pantyhose	●	●	●	●	●	●	●	●
Options									
	Silicone band 3cm	●	●	●	●	●	●	●	●
	Silicone band 5cm	●	●	●	●	●	●	●	●
	Elastic band 5cm	●	●	●	●	●	●	●	●
	Open Toe	●	●	●	●	●	●	●	●
	Closed Toe	●	●	●	●	●	●	●	●
	Oblique G	●	●	●	●	●	●	●	●
	Oblique foot (open or closed toe)	●	●	●	●	●	●	●	●
	Stockings in two parts	●	●	●	●	●	●	●	●
	Shortened instep / T-Heel	●	●	●	●	●	Standard option	Standard option	Standard option
	Full compression body part	-	●	●	●	●	●	●	●
	Curved Knee	●	●	●	●	●	●	●	●
	Horizontal Fly	●	●	●	●	●	●	●	●
	Open Crutch	●	●	●	●	●	●	●	●
Upper Extremity									
AC1	Gauntlet	-	-	-	●	●			
AG	Palm To Axilla	-	-	-	●	●			
CG	Wrist To Axilla	●	●	●	●	●			
Options/Modifications									
	Silicone band 3cm	●	●	●	●	●			
	Silicone band 5cm	●	●	●	●	●			
	Oblique G	●	●	●	-	-			



Support Stockings



Product overview

Medical support stockings.



Return to
contents page

VENOSAN® LEGLINE 20

VENOSAN® LEGLINE 30

ELEGANT LUXURY

- For lighter legs in everyday situations
- Sheer look
- Optimal wearing comfort
- Ensures maximum freedom of movement

Indications




- To prevent varicose veins and heavy legs
- For people who sit or stand a lot at work
- Long journeys by plane, train, and car



Composition
78% polyamide
22% LYCRA® elastane



Composition
AD / AGH:
70% polyamide
30% LYCRA® elastane
AT:
78% polyamide
22% LYCRA® elastane



VENOSAN® MICROFIBERLINE

STYLISH COMFORT, FASHIONABLY SOFT


- Soft and comfortable due to microfibre material
- Stylish support
- Non-restricting top band
- Reinforced heel and toe
- Available in mens and ladies variants

Indications

- Tired, heavy, swollen and aching legs
- Prevention of Deep Vein Thrombosis (DVT)



Composition
Men:
82% Nylon
10% Spandex
8% Cotton
Ladies:
90% Nylon
10% Spandex



VENOSAN® SILVERLINE

EVERYDAY COMFORT

- Enhanced comfort and fit with knitted top band
- Elimination of bacteria and fungus associated with foot odour
- Anti-static effect
- Soft, skin friendly material with ribbed texture
- Wicks away moisture to keep the skin dry and comfortable
- Available in mens and ladies variants

Indications

- Tired, heavy, swollen and aching legs
- Prevention of Deep Vein Thrombosis (DVT)



Composition
Men:
86% Nylon
9% Spandex
5% X-STATIC® Silver plated Nylon
Ladies:
80% Nylon
11% Spandex
9% X-STATIC® Silver plated Nylon



Support Stockings

Product overview

Medical support stockings.



Return to contents page

VENOSAN® ATHLETIC PERFORMANCE SUPPORT

- Unisex compression socks
- Reciprocated heel and toe for improved fit and comfort
- Soft, skin-friendly material with ribbed texture
- Ventilated foot and heel for breathability
- Cushioned achilles tendon area for better protection and comfort

Indications

- Tired, heavy, swollen and aching legs
- Prevention of Deep Vein Thrombosis (DVT)



AD
Below knee stocking

20-30mmHg

Composition
75% Nylon
12% Polyester
13% Spandex



White

Black

Pink

Blue

Purple

VENOSAN® SPORTSLINE PERFORMANCE SUPPORT

- Unisex compression socks
- Made with COOLMAX® fibres that provide superior wicking to keep feet and legs dry
- Cushioned comfort sole with extra padding support
- Wide toe box for roomy comfort
- Easy to put on

Indications

- Tired, heavy, swollen and aching legs
- Prevention of Deep Vein Thrombosis (DVT)



AD
Below knee stocking

15-20mmHg

20-30mmHg

Composition
65% COOLMAX
26% Nylon
9% Spandex



White

Black

Anti-embolism stockings

VENOSAN® AES THROMBOSIS AND EMBOLISM PREVENTION

- Anatomical design providing a good fit
- Colour coding for simple and reliable identification of size
- Top inspection window for safety

Indications

- To prevent thrombosis and embolism before, during and after surgeries
- General immobility, confinement to bed
- Before, during and after childbirth



AD
Below knee stocking



AG
Thigh Length



AGG
Thigh length with belt

18mmHg

Composition
AD / AGG:
91% nylon polyamide
9% LYCRA® elastane
AG:
89% nylon polyamide
11% LYCRA® elastane



White



Return to contents page

VENOSAN® LEGLINE 20

ELEGANT LUXURY



Elegant appearance for your legs, VENOSAN LEGLINE provides effective compression stockings and supports socks with extremely fine mesh.

- For non-heavy legs in everyday situations.
- Retains shape, ensures maximum freedom of movement.
- Sheer look.
- Optimal wearing comfort.
- Approximately 20 mmHg pressure at the ankle.

Lace top band





Outside face



Inside

[Click for VENOSAN LEGLINE 20 ordering and sizing details](#)



20 mm Hg
140 den

Garment Styles



AD
Below knee stocking



AGH
Thigh High Lace top



AT
Pantyhose



ATU
Maternity pantyhose

Indications

- To prevent varicose veins and heavy legs.
- For people who sit or stand a lot at work.
- Long journeys by plane, train, and car.

Composition

78 % Polyamide | 22 % LYCRA® Elastane



Return to contents page

VENOSAN® LEGLINE 30

ELEGANT LUXURY



Elegant appearance for your legs, VENOSAN LEGLINE provides effective compression stockings and supports socks with extremely fine mesh.

- For non-heavy legs in everyday situations.
- Retains shape, ensures maximum freedom of movement.
- Sheer look.
- Optimal wearing comfort.
- Approximately 30 mmHg pressure at the ankle.

Lace top band





Outside face



Inside

[Click for VENOSAN LEGLINE 30 ordering and sizing details](#)



30 mm Hg
280 den

Garment Styles



AD
Below knee stocking



AGH
Thigh High Lace top



AT
Pantyhose



ATU
Maternity pantyhose

Indications

- To prevent varicose veins and heavy legs.
- For people who sit or stand a lot at work.
- Long journeys by plane, train, and car.

Composition

AD | AGH: 70 % Polyamide | 30 % LYCRA® Elastane
AT: 78 % Polyamide | 22 % LYCRA® Elastane



Return to contents page

VENOSAN® MICROFIBERLINE

STYLISH COMFORT, FASHIONABLY SOFT



VENOSAN® MicroFiberLine is a graduated compression support sock with a closed toe. Microfibre softness and excellent breathability ensures wearer comfort. For active men and women in work and travel.


- Graduated compression.
- Stylish comfortable support.
- Fashionably soft comfort.
- Reinforced heel and toe.
- Microfibre softness.
- Comfortable, non-restricting top band.

[Click for VENOSAN MICROFIBERLINE ordering and sizing details](#)

Men



Ladies



Garment Styles



AD
Below knee stocking

Indications

- Tired, heavy, swollen and aching legs.
- Prevention of Deep Vein Thrombosis (DVT).

Composition

Men:
82% Nylon
10% Spandex
8% Cotton

Ladies:
90% Nylon
10% Spandex



Return to contents page

VENOSAN® SILVERLINE

EVERYDAY COMFORT



SILVERLINE® offers long-lasting protection from foot odour and athlete's foot.

X-Static®, the silver fibre, promotes hygiene and anti-microbial properties. Skin remains deodorised and free of bacteria and fungus — particularly effective in protecting against foot odour and athlete's foot. Cool and dry in summer and warm in winter, SILVERLINE® fibres wick moisture away from the skin for comfort in warm weather. X-Static®, the silver fibre, conducts heat in winter to retain warmth. The silver in X-Static® yarns also assists in eliminating static electricity.

Specifically designed

- Enhanced comfort and fit with knitted top band.
- Elimination of bacteria and fungus associated with foot odour.
- Excellent regulator of heat.

[Click for VENOSAN SILVERLINE ordering and sizing details](#)

Men

Black

Ladies

Beige

Black



20–30 mm Hg

Garment Styles



AD
Below knee
stocking

Indications

- Tired, heavy, swollen and aching legs.
- Prevention of Deep Vein Thrombosis (DVT).

Composition

Men:
86% Nylon
9% Spandex
5% X-STATIC®
Silver plated
Nylon

Ladies:
80% Nylon
11% Spandex
9% X-STATIC®
Silver plated
Nylon



Return to contents page

VENOSAN® ATHLETIC PERFORMANCE SUPPORT



Cushioned foot beds and balloon toe pockets, provide day-long benefits of comfort and relief from the effects of venous disease and/or daily living.

- Broad gender appeal.
- Reciprocated heel and toe for outstanding fit and comfort.
- Seamless toe, ventilated foot for breathability and arch support.
- Cushioned achilles tendon area for better protection and comfort.

[Click for VENOSAN ATHLETIC ordering and sizing details](#)



20–30 mm Hg

Garment Styles



AD
Below knee stocking

Indications

- Tired, heavy, swollen and aching legs.
- Prevention of Deep Vein Thrombosis (DVT).

Composition

75% Nylon
12% Polyester
13% Spandex



Return to contents page

VENOSAN® SPORTSLINE

PERFORMANCE SUPPORT




VENOSAN SportsLine is a graduated compression support sock for active men and women. Coolmax® draws moisture away from the skin and helps regulate the body temperature to keep you cool.

Why Coolmax?

Coolmax® is the original performance fabric that has been helping top athletes push harder and longer but is now being discovered for the cool comfort it brings to everyday apparel.

- Cushioned comfort sole with extra padding support.
- Wide toe box or roomy comfort.
- Easy to put on.

[Click for VENOSAN SPORTSLINE ordering and sizing details](#)



Garment Styles



AD
Below knee stocking

Indications

- Tired, heavy, swollen and aching legs.
- Prevention of Deep Vein Thrombosis (DVT).

Composition

65% Coolmax®
26% Nylon
9% Spandex



Lohmann & Rauscher

Anti-Embolic Stockings





Return to contents page

VENOSAN® AES

THROMBOSIS AND EMBOLISM PREVENTION



The VENOSAN AES range comprises a comprehensive selection of well fitting anti thrombosis stockings.

AES stockings work by ensuring constantly decreasing pressure from the ankle region to the knee or thigh, thereby accelerating venous return for patients in the lying position, and helping to prevent thrombosis as a result.

VENOSAN AES stockings are available in three versions: AD, AGH, AGG.

Properties


- Light and air-permeable to ensure maximum comfort.
- Easy to put on.
- Anatomical shaping and comprehensive choice of size.
- Ensuring constantly decreasing pressure from the ankle to the knee or thigh, thereby accelerating venous return in patients in the supine position and helping to prevent thrombosis.
- Top inspection window for safety.
- Can be washed and reused.
- Colour coding at the toes and top seam allows immediate identification of the size.

White

[Click for VENOSAN AES ordering and sizing details](#)

18 mm Hg

Garment Styles



Indications

- To prevent thrombosis and embolism before, during and after surgeries.
- General immobility, confinement to bed.
- Before, during and after childbirth.

Composition

AD / AGH:
91% nylon polyamide
9% LYCRA® elastane

AGG:
89% nylon polyamide
11% LYCRA® elastane



Lohmann & Rauscher

Adjustable Compression Wraps



Click for ReadyWrap ordering and sizing details

 Return to contents page

ReadyWrap™

Adjustable compression garments



ReadyWrap is an adjustable compression wrap that provides patients and their caregivers an easy-to-use solution to support wound healing and oedema management.

ReadyWrap provides clinically effective compression therapy while supporting patient compliance with its unique alternating strapping system.

Product Benefits



Low-stretch material

Mimics short-stretch bandaging — the gold standard of treatment.



50% overlapping straps

Helps prevent gapping.



Blocked support spine

Helps maintain garment form to support all-day wear.



Color-coded straps

Provides easy step-by-step donning guidance.



Foam padding


Cushioning at elbow and wrist provides comfort and evenly distributes pressure.



Modular design

Garments work effectively on their own or in combination.

Colour




Indications

- Lymphoedema
- Chronic oedema
- Venous disease


Available styles



Accessories



Click for
TributeWrap
ordering and
sizing details


Return to
contents page

TributeWrap™

Self-adjustable foam compression garment



TributeWrap is an adjustable compression garment with foam technology that nestles into the tissue to promote skin health.

It supports easy self-management of lymphoedema and oedema during low activity by comfortably adjusting to your patient's unique shape and lifestyle. TributeWrap helps maintain gains made during therapy, while allowing the flexibility to support further reduction.

Simply slide it on, adjust the straps, and cover with Sleep Sleeve for comfortable evening compression.

Product Benefits



Chevron channeling

Designed based on MLD principles of directing fluid to collateral pathways.



Proven foam technology

Breaks up and softens hardened fibrotic tissue to support skin health.



Adjustable strapping system

Easily conforms for the right, comfortable fit.



Easy-to-pull finger grip

Simple pull function to support self-management.



Powered by Coolcore®

Innovative interior fabric with chemical-free Cooling Power which rapidly wicks away sweat and moisture to regulate body temperature.

Colour



Black exterior
Slate interior

Indications

- Mild to severe lymphoedema, stage 1–3
- Chronic oedema (non-blood related)

Available styles






Lohmann & Rauscher

Additional Compression Garments



Click for
Caresia
ordering and
sizing details


Return to
contents page

Caresia™

Off-the-shelf foam bandage liner



Caresia is an off-the-shelf foam bandage liner designed to simplify and enhance compression bandage techniques by replacing the need for traditional padding or under layers. It is to be used during the intensive phase of therapy to simplify the multi-layer bandaging process during the daytime, nighttime, and high or low activity.

Product Benefits



Saves time

Eliminates the underlayers of bandaging. Slide on and apply short-stretch bandaging over top.



Proven foam technology

Breaks up and softens hardened fibrotic tissue to support skin health.



Quick and easy fitting

Only 2–5 measurements needed to select size for each style.



High-stretch fabric

Offers greater range of motion compared to multi-layer bandages.

Colour



Blue

Indications

- Lymphoedema - when used with bandages/other garments
- Oedema (non-blood related) - when used with bandages / other garments

Apply Rosidal K over Caresia to achieve the desired compression profile



Available styles



Gauntlet



MCP to Axilla



Wrist to Axilla



Foot



Glove



Thigh



Below Knee

Click for
ExoSoft
ordering and
sizing details


Return to
contents page

ExoSoft™

Off-the-shelf circular knit compression garments



ExoSoft is an off-the shelf circular-knit gradient compression garment designed for self-management during the day or high activity.

Product Benefits


- Excellent breathability
- Soft and smooth
- Lightweight, breathable microfibre fabric helps keep patients cool, dry and comfortable
- Easy for patients to don
- Wide bands help reduce rolling and keep garment in place
- Armsleeve available with silicone top band
- 6-month warranty
- Easy garment care in washer and dryer

Intended Uses

- Helps maintain reduced volume and shape by supporting extremities
- Management of mild to moderate lymphoedema
- Management of mild to severe oedema (non-blood related)

Colour



Compression	±↑ 20 – 30 mmHg
Top style	KNIT TOP ▒ SILICONE TOP†
Digit style	 QUARTER HALF FULL


† Silicone top is available with Armsleeve

Available styles



Additional Compression Garments

Click for
ExoStrong
ordering and
sizing details


Return to
contents page

ExoStrong™

Off-the-shelf flat knit compression garments



ExoStrong is an off-the-shelf flat-knit gradient compression garment for self-management during the daytime or high activity.

Product Benefits



Supports mobility

High containment garment in a comfortable and functional design.



Quick and easy fitting

Only 2–5 measurements needed to select size for each style.



Off-the-Shelf

Designed based on MLD principles of directing fluid along lymphatic pathways.



Easy garment care

Durable construction ensures pads retain their shape between washes.



Made in the USA

High-quality flat knit manufactured at our headquarters in Milwaukee, WI.

Colour



Beige

Black

Indications

- Mild to severe lymphoedema
- Oedema

Compression	↕ 20–30 mmHg
Top style	KNIT TOP* ▨ SILICONE TOP†
Toe style	○ OPEN TOE
Digit style	☞ QUARTER HALF FULL

* Knit top is available with Arm and Below Knee

† Silicone top is available with Arm, Below Knee, and Thigh High

Available styles




Click for
SwellSpots
ordering and
sizing details


Return to
contents page

SwellSpots™

Oedema and fibrosis pads



SwellSpots are reusable therapeutic foam pads to incorporate into multi-layer bandaging or garments to provide focal compression to help soften troublesome areas of oedema and fibrosis.

Ideal for difficult areas of oedema, fibrosis, skin folds, misshapened limbs as well as padding and protection.

Product Benefits



Proven foam technology

Breaks up and softens hardened fibrotic tissue to support skin health.



Versatile shapes and sizes

Wide range of off-the-shelf and custom solutions to manage troublesome areas no matter the location.



Channeling

Designed based on MLD principles of directing fluid along lymphatic pathways.



Reusable pads

Durable construction ensures pads retain their shape between washes.



Easy care

Machine washable and dryer.

Colour

Slate

Available off-the-shelf styles



Accessories for compression stockings





Return to contents page

Accessories for compression stockings.



VENO-GLIDER

Foot Slip Donning Aid

Article no.	Description	Unit
995345	VENO-GLIDER	1 pc.



VENOSAN® Dot Gloves

Washable gloves for donning and doffing compression stockings, available in two sizes.

Material composition

Knitted fabric: 100% nylon polyamide, Dots: PVC

Article no.	Description	Unit
VC003	VENOSAN Dot Gloves Size S/M	1 pair
VC004	VENOSAN Dot Gloves Size L/XL	1 pair



COLTEX

Special adhesive for securing medical compression stockings to the skin.

Composition

PVP/VA copolymer, glycerin, 2-propanol

Article no.	Description	Unit
KC67000	COLTEX Adhesive 60ml	1 pc.



VENOSAN® Eureka On

Donning aid for compression stockings. Can be used with open and closed toe.

Article no.	Description	Unit
155380	VENOSAN Eureka On	1 pc.



Lohmann & Rauscher

Measuring and ordering information






Return to contents page

VENOSAN® Measuring Guide






When? Where possible, it is recommended to measure in the morning, when legs are rested and swelling is at a minimum.

How? Measurements to be taken using a tape measure (cm only) with the patient standing. For Below knee garments, the patient can be measured when sitting. It is recommended that measurements are taken on a bare leg.

Measuring points



Which measurements are needed

 <p>AD Below knee</p> <p>ADH Below knee with silicone</p> <ul style="list-style-type: none"> ■ cB ■ cC ■ cD ■ lD 	 <p>AGH Thigh length with silicone</p> <ul style="list-style-type: none"> ■ cB ■ cC ■ cD ■ cG ■ lG 	 <p>AGG Thigh length with belt</p> <ul style="list-style-type: none"> ■ cB ■ cC ■ cD ■ cG ■ lK 	 <p>AT Panty hose</p> <ul style="list-style-type: none"> ■ cB ■ cC ■ cD ■ cG ■ lK 	 <p>ATU Maternity pantyhose</p> <ul style="list-style-type: none"> ■ cB ■ cC ■ cD ■ cG ■ lK
---	---	---	---	--

VENOSAN® Sizing Guide

VENOSAN®

Compression Stockings



Return to
contents page

VENOSAN 3000 **VENOSAN 4000** **VENOSAN 5000** **VENOSAN 7000** **VENOSAN 8000**

Circumference

	Small (cm)	Medium (cm)	Large (cm)	X-Large (cm)
cB - Ankle	20 – 23	23 – 26	26 – 29	29 – 31
cC - Calf	29 – 39	33 – 42	36 – 45	39 – 47
cD - Below knee	28 – 38	32 – 41	35 – 44	38 – 46
cG - Thigh	48 – 63	54 – 67	58 – 72	63 – 74

Length

	Length of stocking	Short (cm)	Long (cm)
ℓD - Below knee	AD / ADH	35 – 40	39 – 44
ℓG - Below gluteal fold	AGH	63 – 72	68 – 78
ℓK - Above gluteal fold	AGG / AT / ATU Maternity	67 – 74	74 – 85

VENOSAN 6000

Circumference

	X-Small (cm)	Small (cm)	Medium (cm)	Large (cm)	X-Large (cm)
cB - Ankle	16 – 19	20 – 22	23 – 25	26 – 28	29 – 31
cC - Calf	27 – 36	30 – 39	33 – 41	36 – 44	39 – 47
cD - Below knee	26 – 35	29 – 38	32 – 40	35 – 43	38 – 46
cG - Thigh	44 – 59	48 – 63	54 – 66	58 – 70	63 – 74

Length

	Length of stocking	Short (cm)	Long (cm)
ℓD - Below knee	AD / ADH	34 – 39	39 – 44
ℓG - Below gluteal fold	AGH	59 – 67	68 – 78
ℓK - Above gluteal fold	AGG / AT	67 – 74	74 – 85

If patient does not fit into the size range provided, we also have Made To Measure (MTM) available. Contact L&R on 1300 572 869 for more information.





Return to contents page


VENOSAN®

Measuring and Sizing Guide for Arm Compression Stockings

When? It is recommended to measure the arm when it is in a relaxed position. The length is measured on the outstretched arm.

How? Measurements to be taken using a tape measure (cm only).

Measuring Points



Sizing Guide

Circumference wrist

	Small (cm)	Medium (cm)	Large (cm)	X-Large (cm)
cC	14 – 17	17 – 20	20 – 23	23 – 26

Circumference forearm

	Small (cm)	Medium (cm)	Large (cm)	X-Large (cm)
cD	22 – 26	26 – 30	30 – 34	34 – 38


Circumference upper arm

	Small (cm)	Medium (cm)	Large (cm)	X-Large (cm)
cF	24 – 30	30 – 36	36 – 42	42 – 48


Length

	Short (cm)	Long (cm)
Wrist (C) to Axilla (G)	35 – 45	40 – 50

Styles available



CG
Wrist to Axilla



CG-H
Wrist to Axilla with silicone band

VENOSAN® 3000 COTTON NATURAL LOOK

Click for
VENOSAN 3000
product page



Return to
contents page

off-the-shelf range

VENOSAN® 3000 COTTON NATURAL LOOK

CCL
2

Style	Colour	Length	Toe	Silicone Band	REF			
					S	M	L	XL
AD Below knee	beige	long	open	–	SB32001	SB32002	SB32003	SB32004
		long	closed	–	SB32011	SB32012	SB32013	SB32014
		short	open	–	SB32021	SB32022	SB32023	SB32024
		short	closed	–	SB32031	SB32032	SB32033	SB32034
ADH Below knee with silicone	beige	long	open	plain top	SB32041	SB32042	SB32043	SB32044
		long	closed	plain top	SB32051	SB32052	SB32053	SB32054
		short	open	plain top	SB32061	SB32062	SB32063	SB32064
		short	closed	plain top	SB32071	SB32072	SB32073	SB32074
AGH Thigh length with silicone	beige	long	open	plain top	SB32601	SB32602	SB32603	SB32604
		short	open	plain top	SB32621	SB32622	SB32623	SB32624

VENOSAN® 3000 COTTON NATURAL LOOK

CCL
3

Style	Colour	Length	Toe	Silicone Band	REF			
					S	M	L	XL
AD Below knee	beige	long	open	–	SB33001	SB33002	SB33003	SB33004
		short	open	–	SB33021	SB33022	SB33023	SB33024

VENOSAN® 4000

CLIMATE EFFECT

Click for
VENOSAN 4000
product page



Return to
contents page

off-the-shelf range

VENOSAN® 4000 CLIMATE EFFECT

CCL
1

Style	Colour	Length	Toe	Silicone Band	REF			
					S	M	L	XL
AD Below knee	mexico	long	open	–	SG41001	SG41002	SG41003	SG41004
		long	closed	–	SG41011	SG41012	SG41013	SG41014
		short	open	–	SG41021	SG41022	SG41023	SG41024
		short	closed	–	SG41031	SG41032	SG41033	SG41034
	marokko	long	open	–	SM41001	SM41002	SM41003	SM41004
		long	closed	–	SM41011	SM41012	SM41013	SM41014
		short	open	–	SM41021	SM41022	SM41023	SM41024
		short	closed	–	SM41031	SM41032	SM41033	SM41034
	black	long	open	–	SN41001	SN41002	SN41003	SN41004
		long	closed	–	SN41011	SN41012	SN41013	SN41014
		short	open	–	SN41021	SN41022	SN41023	SN41024
		short	closed	–	SN41031	SN41032	SN41033	SN41034
ADH Below knee with silicone	mexico	long	open	plain top	SG41041	SG41042	SG41043	SG41044
		long	closed	plain top	SG41051	SG41052	SG41053	SG41054
		short	open	plain top	SG41061	SG41062	SG41063	SG41064
		short	closed	plain top	SG41071	SG41072	SG41073	SG41074
	marokko	long	open	plain top	SM41041	SM41042	SM41043	SM41044
		long	closed	plain top	SM41051	SM41052	SM41053	SM41054
		short	open	plain top	SM41061	SM41062	SM41063	SM41064
		short	closed	plain top	SM41071	SM41072	SM41073	SM41074
	black	long	open	plain top	SN41041	SN41042	SN41043	SN41044
		long	closed	plain top	SN41051	SN41052	SN41053	SN41054
		short	open	plain top	SN41061	SN41062	SN41063	SN41064
		short	closed	plain top	SN41071	SN41072	SN41073	SN41074
AGH Thigh length with silicone	mexico	long	closed	lace top	SG41601	SG41602	SG41603	SG41604
		long	open	plain top	SG41611	SG41612	SG41613	SG41614
		long	closed	plain top	SG41621	SG41622	SG41623	SG41624
		short	closed	lace top	SG41631	SG41632	SG41633	SG41634
		short	open	plain top	SG41641	SG41642	SG41643	SG41644
		short	closed	plain top	SG41651	SG41652	SG41653	SG41654
	marokko	long	closed	lace top	SM41601	SM41602	SM41603	SM41604
		long	open	plain top	SM41611	SM41612	SM41613	SM41614
		long	closed	plain top	SM41621	SM41622	SM41623	SM41624
		short	closed	lace top	SM41631	SM41632	SM41633	SM41634
		short	open	plain top	SM41641	SM41642	SM41643	SM41644
		short	closed	plain top	SM41651	SM41652	SM41653	SM41654
black	long	closed	lace top	SN41601	SN41602	SN41603	SN41604	
	long	open	plain top	SN41611	SN41612	SN41613	SN41614	
	long	closed	plain top	SN41621	SN41622	SN41623	SN41624	
	short	closed	lace top	SN41631	SN41632	SN41633	SN41634	
	short	open	plain top	SN41641	SN41642	SN41643	SN41644	
	short	closed	plain top	SN41651	SN41652	SN41653	SN41654	

VENOSAN® 4000 CLIMATE EFFECT

Click for
VENOSAN 4000
product page



Return to
contents page

off-the-shelf range

VENOSAN® 4000 CLIMATE EFFECT

CCL
1

Style	Colour	Length	Toe	Silicone Band	REF			
					S	M	L	XL
AT Pantyhose	mexico	long	closed	–	SG41221	SG41222	SG41223	SG41224
		short	closed	–	SG41231	SG41232	SG41233	SG41234
	marokko	long	closed	–	SM41221	SM41222	SM41223	SM41224
		short	closed	–	SM41231	SM41232	SM41233	SM41234
	black	long	closed	–	SN41221	SN41222	SN41223	SN41224
		short	closed	–	SN41231	SN41232	SN41233	SN41234
ATU Maternity Pantyhose	mexico	long	closed	–	SG41301	SG41302	SG41303	SG41304
	black	long	closed	–	SN41301	SN41302	SN41303	SN41304

VENOSAN® 4000 CLIMATE EFFECT

CCL
2

Style	Colour	Length	Toe	Silicone Band	REF			
					S	M	L	XL
AD Below knee	mexico	long	open	–	SG42001	SG42002	SG42003	SG42004
		long	closed	–	SG42011	SG42012	SG42013	SG42014
		short	open	–	SG42021	SG42022	SG42023	SG42024
		short	closed	–	SG42031	SG42032	SG42033	SG42034
	marokko	long	open	–	SM42001	SM42002	SM42003	SM42004
		long	closed	–	SM42011	SM42012	SM42013	SM42014
		short	open	–	SM42021	SM42022	SM42023	SM42024
		short	closed	–	SM42031	SM42032	SM42033	SM42034
	black	long	open	–	SN42001	SN42002	SN42003	SN42004
		long	closed	–	SN42011	SN42012	SN42013	SN42014
		short	open	–	SN42021	SN42022	SN42023	SN42024
		short	closed	–	SN42031	SN42032	SN42033	SN42034
ADH Below knee with silicone	mexico	long	open	plain top	SG42041	SG42042	SG42043	SG42044
		long	closed	plain top	SG42051	SG42052	SG42053	SG42054
		short	open	plain top	SG42061	SG42062	SG42063	SG42064
		short	closed	plain top	SG42071	SG42072	SG42073	SG42074
	marokko	long	open	plain top	SM42041	SM42042	SM42043	SM42044
		long	closed	plain top	SM42051	SM42052	SM42053	SM42054
		short	open	plain top	SM42061	SM42062	SM42063	SM42064
		short	closed	plain top	SM42071	SM42072	SM42073	SM42074
	black	long	open	plain top	SN42041	SN42042	SN42043	SN42044
		long	closed	plain top	SN42051	SN42052	SN42053	SN42054
		short	open	plain top	SN42061	SN42062	SN42063	SN42064
		short	closed	plain top	SN42071	SN42072	SN42073	SN42074

VENOSAN® 4000 CLIMATE EFFECT

Click for
VENOSAN 4000
product page



Return to
contents page

off-the-shelf range

CCL
2

Style	Colour	Length	Toe	Silicone Band	REF			
					S	M	L	XL
AGH Thigh length with silicone	mexico	long	closed	lace top	SG42601	SG42602	SG42603	SG42604
		long	open	plain top	SG42611	SG42612	SG42613	SG42614
		long	closed	plain top	SG42621	SG42622	SG42623	SG42624
		short	closed	lace top	SG42631	SG42632	SG42633	SG42634
		short	open	plain top	SG42641	SG42642	SG42643	SG42644
		short	closed	plain top	SG42651	SG42652	SG42653	SG42654
	marokko	long	closed	lace top	SM42601	SM42602	SM42603	SM42604
		long	open	plain top	SM42611	SM42612	SM42613	SM42614
		long	closed	plain top	SM42621	SM42622	SM42623	SM42624
		short	closed	lace top	SM42631	SM42632	SM42633	SM42634
		short	open	plain top	SM42641	SM42642	SM42643	SM42644
		short	closed	plain top	SM42651	SM42652	SM42653	SM42654
	black	long	closed	lace top	SN42601	SN42602	SN42603	SN42604
		long	open	plain top	SN42611	SN42612	SN42613	SN42614
		long	closed	plain top	SN42621	SN42622	SN42623	SN42624
		short	closed	lace top	SN42631	SN42632	SN42633	SN42634
		short	open	plain top	SN42641	SN42642	SN42643	SN42644
		short	closed	plain top	SN42651	SN42652	SN42653	SN42654
AGG (left) Thigh length with belt	mexico	long	open	-	SG42151	SG42152	SG42153	SG42154
AGG (right) Thigh length with belt	mexico	long	open	-	SG42171	SG42172	SG42173	SG42174
AT Pantyhose	mexico	long	closed	-	SG42221	SG42222	SG42223	SG42224
		short	closed	-	SG42231	SG42232	SG42233	SG42234
	marokko	long	closed	-	SM42221	SM42222	SM42223	SM42224
		short	closed	-	SM42231	SM42232	SM42233	SM42234
	black	long	closed	-	SN42221	SN42222	SN42223	SN42224
		short	closed	-	SN42231	SN42232	SN42233	SN42234
ATU Maternity Pantyhose	mexico	long	closed	-	SG42301	SG42302	SG42303	SG42304
		short	closed	-	SG42311	SG42312	SG42313	SG42314
	black	long	closed	-	SN42301	SN42302	SN42303	SN42304
		long	closed	-	SN42311	SN42312	SN42313	SN42314
		short	closed	-	SN42301	SN42302	SN42303	SN42304
		short	closed	-	SN42311	SN42312	SN42313	SN42314

made-to-measure

VENOSAN® 4000 CLIMATE EFFECT

CCL
1

CCL
2

Style	Colour	Toe	CCL 1	CCL 2
AD	CCL 1 CCL 2 mexico, black	open/closed	•	•
ADH		open/closed	•	•
AG		open/closed	•	•
AGH		open/closed	•	•
AGG		open/closed	•	•
AT		open/closed	•	•
ATU Maternity		open/closed	•	•
		open/closed	•	•

• = available

VENOSAN® 5000

THE MAGIC OF SILVER

Click for
VENOSAN 5000
product page



Return to
contents page

off-the-shelf range

VENOSAN® 5000

THE MAGIC OF SILVER

CCL
1

Style	Colour	Length	Toe	Silicone Band	REF			
					S	M	L	XL
AD Below knee	mexico	long	open	–	SG45001	SG45002	SG45003	SG45004
		long	closed	–	SG45011	SG45012	SG45013	SG45014
	marokko	long	open	–	SM45001	SM45002	SM45003	SM45004
		long	closed	–	SM45011	SM45012	SM45013	SM45014
	black	long	open	–	SN45001	SN45002	SN45003	SN45004
long		closed	–	SN45011	SN45012	SN45013	SN45014	
AGH Thigh length with silicone	mexico	long	closed	lace top	SG45601	SG45602	SG45603	SG45604
		long	open	plain top	SG45611	SG45612	SG45613	SG45614
	marokko	long	closed	lace top	SM45601	SM45602	SM45603	SM45604
		long	open	plain top	SM45611	SM45612	SM45613	SM45614
	black	long	closed	lace top	SN45601	SN45602	SN45603	SN45604
		long	open	plain top	SN45611	SN45612	SN45613	SN45614
AGG (left) Thigh length with belt	mexico	long	open	–	SG45151	SG45152	SG45153	SG45154
AGG (right) Thigh length with belt	mexico	long	open	–	SG45171	SG45172	SG45173	SG45174
AT Pantyhose	mexico	long	closed	–	SG45221	SG45222	SG45223	SG45224
	marokko	long	closed	–	SM45221	SM45222	SM45223	SM45224
	black	long	closed	–	SN45221	SN45222	SN45223	SN45224

VENOSAN® 5000

THE MAGIC OF SILVER


CCL
2

Style	Colour	Length	Toe	Silicone Band	REF			
					S	M	L	XL
AD Below knee	mexico	long	open	–	SG46001	SG46002	SG46003	SG46004
		long	closed	–	SG46011	SG46012	SG46013	SG46014
	marokko	long	open	–	SM46001	SM46002	SM46003	SM46004
		long	closed	–	SM46011	SM46012	SM46013	SM46014
	black	long	open	–	SN46001	SN46002	SN46003	SN46004
long		closed	–	SN46011	SN46012	SN46013	SN46014	
AGH Thigh length with silicone	mexico	long	open	plain top	SG46611	SG46612	SG46613	SG46614
		long	closed	plain top	SG46621	SG46622	SG46623	SG46624
	marokko	long	open	plain top	SM46611	SM46612	SM46613	SM46614
		long	closed	plain top	SM46621	SM46622	SM46623	SM46624
	black	long	open	plain top	SN46611	SN46612	SN46613	SN46614
		long	closed	plain top	SN46621	SN46622	SN46623	SN46624
AGG (left) Thigh length with belt	mexico	long	open	–	SG46151	SG46152	SG46153	SG46154
AGG (right) Thigh length with belt	mexico	long	open	–	SG46171	SG46172	SG46173	SG46174
AT Pantyhose	mexico	long	closed	–	SG46221	SG46222	SG46223	SG46224
	marokko	long	closed	–	SM46221	SM46222	SM46223	SM46224
	black	long	closed	–	SN46221	SN46222	SN46223	SN46224

VENOSAN® 5000

THE MAGIC OF SILVER

Click for
VENOSAN 5000
product page



off-the-shelf range

VENOSAN® 5000 THE MAGIC OF SILVER

CCL
3

Style	Colour	Length	Toe	Silicone Band	REF			
					S	M	L	XL
AD Below knee	mexico	long	open	–	SG47001	SG47002	SG47003	SG47004
		long	closed	–	SG47011	SG47012	SG47013	SG47014
	marokko	long	open	–	SM47001	SM47002	SM47003	SM47004
		long	closed	–	SM47011	SM47012	SM47013	SM47014
	black	long	open	–	SN47001	SN47002	SN47003	SN47004
		long	closed	–	SN47011	SN47012	SN47013	SN47014

made-to-measure

VENOSAN® 5000 THE MAGIC OF SILVER

CCL
1

CCL
2

CCL
3

Style	Colour	Toe	CCL 1	CCL 2	CCL 3
AD	CCL 1 CCL 2 mexico, black CCL 3 mexico	open/closed	•	•	•
ADH		open/closed	•	•	•
AG		open/closed	•	•	•
AGH		open/closed	•	•	•
AGG		open/closed	•	•	•
AT		open/closed	•	•	•
ATU Maternity		open/closed	•	•	•

• = available

VENOSAN® 6000 CLASSIC

Click for
VENOSAN 6000
product page



Return to
contents page

off-the-shelf range

VENOSAN® 6000 CLASSIC

CCL
1

Style	Colour	Length	Toe	Silicone Band	REF				
					XS	S	M	L	XL
AD Below knee	beige	long	open	–	CG61000	CG61001	CG61002	CG61003	CG61004
		long	closed	–	CG61010	CG61011	CG61012	CG61013	CG61014
		short	open	–	CG61020	CG61021	CG61022	CG61023	CG61024
		short	closed	–	CG61030	CG61031	CG61032	CG61033	CG61034
	black	long	open	–	CN61000	CN61001	CN61002	CN61003	CN61004
		long	closed	–	CN61010	CN61011	CN61012	CN61013	CN61014
		short	open	–	CN61020	CN61021	CN61022	CN61023	CN61024
		short	closed	–	CN61030	CN61031	CN61032	CN61033	CN61034
AGH Thigh length with silicone	beige	long	open	plain top	CG61600	CG61601	CG61602	CG61603	CG61604
		long	closed	plain top	CG61610	CG61611	CG61612	CG61613	CG61614
		short	open	plain top	CG61620	CG61621	CG61622	CG61623	CG61624
		short	closed	plain top	CG61630	CG61631	CG61632	CG61633	CG61634
	black	long	open	plain top	CN61600	CN61601	CN61602	CN61603	CN61604
		long	closed	plain top	CN61610	CN61611	CN61612	CN61613	CN61614
		short	open	plain top	CN61620	CN61621	CN61622	CN61623	CN61624
		short	closed	plain top	CN61630	CN61631	CN61632	CN61633	CN61634
AT Pantyhose	beige	long	closed	–	CG61210	CG61211	CG61212	CG61213	CG61214
		short	closed	–	CG61230	CG61231	CG61232	CG61233	CG61234
	black	long	closed	–	CN61210	CN61211	CN61212	CN61213	CN61214
		short	closed	–	CN61230	CN61231	CN61232	CN61233	CN61234

VENOSAN® 6000 CLASSIC

CCL
2

Style	Colour	Length	Toe	Silicone Band	REF				
					XS	S	M	L	XL
AD Below knee	beige	long	open	–	CG62000	CG62001	CG62002	CG62003	CG62004
		long	closed	–	CG62010	CG62011	CG62012	CG62013	CG62014
		short	open	–	CG62020	CG62021	CG62022	CG62023	CG62024
		short	closed	–	CG62030	CG62031	CG62032	CG62033	CG62034
	black	long	open	–	CN62000	CN62001	CN62002	CN62003	CN62004
		long	closed	–	CN62010	CN62011	CN62012	CN62013	CN62014
		short	open	–	CN62020	CN62021	CN62022	CN62023	CN62024
		short	closed	–	–	CN62231	CN62232	CN62233	CN62234

VENOSAN® 6000 CLASSIC

Click for
VENOSAN 6000
product page



Return to
contents page

off-the-shelf range

VENOSAN® 6000 CLASSIC

CCL
2

Style	Colour	Length	Toe	Silicone Band	REF				
					XS	S	M	L	XL
AGH Thigh length with silicone	beige	long	open	plain top	CG62600	CG62601	CG62602	CG62603	CG62604
		long	closed	plain top	–	CG62611	CG62612	CG62613	CG62614
		short	open	plain top	CG62620	CG62621	CG62622	CG62623	CG62624
		short	closed	plain top	–	CG62631	CG62632	CG62633	CG62634
	black	long	open	plain top	CN62600	CN62601	CN62602	CN62603	CN62604
		long	closed	plain top	–	CN62611	CN62612	CN62613	CN62614
		short	open	plain top	CN62620	CN62621	CN62622	CN62623	CN62624
		short	closed	plain top	–	CN62631	CN62632	CN62633	CN62634
AGG (left) Thigh length with belt	beige	long	open	–	CG62150	CG62151	CG62152	CG62153	CG62154
		short	open	–	CG62160	CG62161	CG62162	CG62163	CG62164
AGG (right) Thigh length with belt	beige	long	open	–	CG62170	CG62171	CG62172	CG62173	CG62174
		short	open	–	CG62180	CG62181	CG62182	CG62183	CG62184
AT Pantyhose	beige	long	closed	–	CG62210	CG62211	CG62212	CG62213	CG62214
		short	closed	–	CG62230	CG62231	CG62232	CG62233	CG62234
	black	long	closed	–	CN62210	CN62211	CN62212	CN62213	CN62214
		short	closed	–	CN62230	CN62231	CN62232	CN62233	CN62234

VENOSAN® 6000 CLASSIC

CCL
3

Style	Colour	Length	Toe	Silicone Band	REF				
					XS	S	M	L	XL
AD Below knee	beige	long	open	–	CG63000	CG63001	CG63002	CG63003	CG63004
		short	open	–	CG63020	CG63021	CG63022	CG63023	CG63024
AGH Thigh length with silicone	beige	long	open	plain top	CG63600	CG63601	CG63602	CG63603	CG63604
		short	open	plain top	CG63620	CG63621	CG63622	CG63623	CG63624
AT Pantyhose	beige	long	open	–	CG63200	CG63201	CG63202	CG63203	CG63204
		short	open	–	CG63020	CG63021	CG63022	CG63023	CG63024

VENOSAN® 7000 NATURAL TOUCH

Click for
VENOSAN 7000
product page



Return to
contents page

off-the-shelf range

VENOSAN® 7000 NATURAL TOUCH

CCL
1

Style	Colour	Length	Toe	Silicone Band	REF			
					S	M	L	XL
AD Below knee	beige	long	open	–	SG71001	SG71002	SG71003	SG71004
		long	closed	–	SG71011	SG71012	SG71013	SG71014
		short	open	–	SG71021	SG71022	SG71023	SG71024
		short	closed	–	SG71031	SG71032	SG71033	SG71034

VENOSAN® 7000 NATURAL TOUCH

CCL
2

Style	Colour	Length	Toe	Silicone Band	REF			
					S	M	L	XL
AD Below knee	beige	long	open	–	SG72001	SG72002	SG72003	SG72004
		long	closed	–	SG72011	SG72012	SG72013	SG72014
		short	open	–	SG72021	SG72022	SG72023	SG72024
		short	closed	–	SG72031	SG72032	SG72033	SG72034
	black	long	open	–	SN72001	SN72002	SN72003	SN72004
		long	closed	–	SN72011	SN72012	SN72013	SN72014
		short	open	–	SN72021	SN72022	SN72023	SN72024
		short	closed	–	SN72031	SN72032	SN72033	SN72034
ADH Below knee with silicone	beige	long	open	plain top	SG72041	SG72042	SG72043	SG72044
		long	closed	plain top	SG72051	SG72052	SG72053	SG72054
		short	open	plain top	SG72061	SG72062	SG72063	SG72064
		short	closed	plain top	SG72071	SG72072	SG72073	SG72074
	black	long	open	plain top	SN72041	SN72042	SN72043	SN72044
		long	closed	plain top	SN72051	SN72052	SN72053	SN72054
		short	open	plain top	SN72061	SN72062	SN72063	SN72064
		short	closed	plain top	SN72071	SN72072	SN72073	SN72074
AGH Thigh length with silicone	beige	long	open	plain top	SG72601	SG72602	SG72603	SG72604
		long	closed	plain top	SG72611	SG72612	SG72613	SG72614
		short	open	plain top	SG72621	SG72622	SG72623	SG72624
		short	closed	plain top	SG72631	SG72632	SG72633	SG72634
	black	long	open	plain top	SN72601	SN72602	SN72603	SN72604
		long	closed	plain top	SN72611	SN72612	SN72613	SN72614
		short	open	plain top	SN72621	SN72622	SN72623	SN72624
		short	closed	plain top	SN72631	SN72632	SN72633	SN72634
AGG (left) Thigh length with belt	beige	long	open	–	SG72151	SG72152	SG72153	SG72154
		short	open	–	SG72161	SG72162	SG72163	SG72164
AGG (right) Thigh length with belt	beige	long	open	–	SG72171	SG72172	SG72173	SG72174
		short	open	–	SG72181	SG72182	SG72183	SG72184

VENOSAN® 7000 NATURAL TOUCH

Click for
VENOSAN 7000
product page



Return to
contents page

off-the-shelf range

VENOSAN® 7000 NATURAL TOUCH

CCL
3

Style	Colour	Length	Toe	Silicone Band	REF			
					S	M	L	XL
AD Below knee	beige	long	open	–	SG73001	SG73002	SG73003	SG73004
		long	closed	–	SG73011	SG73012	SG73013	SG73014
		short	open	–	SG73021	SG73022	SG73023	SG73024
		short	closed	–	SG73031	SG73032	SG73033	SG73034
	black	long	open	–	SN73001	SN73002	SN73003	SN73004
		long	closed	–	SN73011	SN73012	SN73013	SN73014
		short	open	–	SN73021	SN73022	SN73023	SN73024
		short	closed	–	SN73031	SN73032	SN73033	SN73034
ADH Below knee with silicone	beige	long	open	plain top	SG73041	SG73042	SG73043	SG73044
		long	closed	plain top	SG73051	SG73052	SG73053	SG73054
		short	open	plain top	SG73061	SG73062	SG73063	SG73064
		short	closed	plain top	SG73071	SG73072	SG73073	SG73074
	black	long	open	plain top	SN73041	SN73042	SN73043	SN73044
		long	closed	plain top	SN73051	SN73052	SN73053	SN73054
		short	open	plain top	SN73061	SN73062	SN73063	SN73064
		short	closed	plain top	SN73071	SN73072	SN73073	SN73074
AGH Thigh length with silicone	beige	long	open	plain top	SG73601	SG73602	SG73603	SG73604
		long	closed	plain top	SG73611	SG73612	SG73613	SG73614
		short	open	plain top	SG73621	SG73622	SG73623	SG73624
		short	closed	plain top	SG73631	SG73632	SG73633	SG73634
	black	long	open	plain top	SN73601	SN73602	SN73603	SN73604
		long	closed	plain top	SN73611	SN73612	SN73613	SN73614
		short	open	plain top	SN73621	SN73622	SN73623	SN73624
		short	closed	plain top	SN73631	SN73632	SN73633	SN73634

made-to-measure

VENOSAN® 7000 NATURAL TOUCH

CCL
2

CCL
3

Style	Colour	Toe	CCL 2	CCL 3
AD	beige, black	open/closed	•	•
ADH		open/closed	•	•
AG		open/closed	•	•
AGH		open/closed	•	•
AGG		open/closed	•	•
AT		open/closed	•	•
ATH		open/closed	•	•

• = available

VENOSAN® 8000 ULCERFIT COMPRESSION STOCKING SYSTEM

Click for
VENOSAN 8000
product page



Return to
contents page

off-the-shelf range

VENOSAN® 8000 ULCERFIT COMPRESSION STOCKING SYSTEM

CCL
2

Style	Colour	Length	Toe	Silicone Band	REF			
					S	M	L	XL
AD (Set*) Below knee	beige / white	long 37 – 42cm	understocking: closed	–	SG82001	SG82002	SG82003	SG82004
		short 35 – 40cm	overstocking: open	–	SG82021	SG82022	SG82023	SG82024

* Set: consisting of one (beige) overstocking and two (white) understockings

VENOSAN® 8000 ULCERFIT COMPRESSION STOCKING SYSTEM

CCL
3

Style	Colour	Length	Toe	Silicone Band	REF			
					S	M	L	XL
AD (Set*) Below knee	beige / white	long 37 – 42cm	understocking: closed	–	SG83001	SG83002	SG83003	SG83004
		short 35 – 40cm	overstocking: open	–	SG83021	SG83022	SG83023	SG83024

* Set: consisting of one (beige) overstocking and two (white) understockings

VENOSAN® 8000 ULCERFIT LINERS

CCL
3

Style	Colour	Length	Toe	Silicone Band	REF			
					S	M	L	XL
AD (Pair) Below knee	white	35 – 42cm	closed	–	SG80001	SG80002	SG80003	SG80004

VENOSAN® 4000 Armsleeves

CLIMATE EFFECT

Click for
VENOSAN 4000
Armsleeves
product page



Return to
contents page

off-the-shelf range

VENOSAN® 4000 Armsleeves CLIMATE EFFECT

CCL
2

Style	Colour	Length	Silicone Band	REF			
				S	M	L	XL
CG Wrist to axilla	mexico	long	–	SG42501	SG42502	SG42503	SG42504

made-to-measure

VENOSAN® 4000 Armsleeves CLIMATE EFFECT

Style	Colour	With or without hand piece	4002 CCL 2
CG	beige, black	without hand	•
CG-H		without hand	•

• = available

VENOSAN® 7000 Armsleeves

NATURAL TOUCH

Click for
VENOSAN 7000
Armsleeves
product page



Return to
contents page

off-the-shelf range

VENOSAN® 7000 Armsleeves

NATURAL TOUCH

CCL
2

Style	Colour	Length	With or without hand piece	Silicone Band	REF			
					S	M	L	XL
CG Wrist to axilla	beige	long	without hand	–	SG72501	SG72502	SG72503	SG72504
		short	without hand	–	SG72511	SG72512	SG72513	SG72514
CG-H Wrist to axilla with silicone	beige	long	without hand	plain top	SG72521	SG72522	SG72523	SG72524
		short	without hand	plain top	SG72531	SG72532	SG72533	SG72534
ZG MCP to axilla	beige	long	with hand	–	SG72541	SG72542	SG72543	SG72544
		short	with hand	–	SG72551	SG72552	SG72553	SG72554
ZG-H MCP to Axilla with silicone	beige	long	with hand	plain top	SG72561	SG72562	SG72563	SG72564
		short	with hand	plain top	SG72571	SG72572	SG72573	SG72574

made-to-measure

VENOSAN® 7000 Armsleeves

NATURAL TOUCH

Style	Colour	With or without hand piece	7002 CCL 2	7003 CCL 3
CG	beige, black	without hand	•	•
CG-H		without hand	•	•
CH		without hand	•	•
ZG		with hand	•	•
ZG-H		with hand	•	•
ZH		with hand	•	•

• = available

VENOSAN® LEGLINE 20

ELEGANT LUXURY

Click for
VENOSAN
LEGLINE 20
product page




Return to
contents page

off-the-shelf range

Sizing Guide

Circumference	X-Small (cm)	Small (cm)	Medium (cm)	Large (cm)	X-Large (cm)
cB - Ankle	20 – 21	22 – 23	24 – 25	26 – 27	28 – 29
cC - Calf	31 – 33	34 – 36	37 – 39	40 – 42	43 – 45
cD - Below knee	29 – 31	32 – 34	35 – 37	38 – 40	41 – 43
cF - Mid thigh	42 – 44	45 – 47	48 – 52	51 – 53	54 – 60



Ordering Information

VENOSAN® LEGLINE 20

ELEGANT LUXURY

20mmHg

Style	Colour	Toe	Silicone Band	REF				
				XS	S	M	L	XL
AD Below knee	nude	closed	–	VL303NU	VL304NU	VL305NU	VL306NU	
	sahara	closed	–	VL303SA	VL304SA	VL305SA	VL306SA	
	black	closed	–	VL303BL	VL304BL	VL305BL	VL306BL	
AGH Thigh length with silicone	nude	closed	lace top	VL323NU	VL324NU	VL325NU	VL326NU	
	sahara	closed	lace top	VL323SA	VL324SA	VL325SA	VL326SA	
	black	closed	lace top	VL323BL	VL324BL	VL325BL	VL326BL	
AT Pantyhose	nude	closed	–	VL333NU	VL334NU	VL335NU	VL336NU	VL337NU
	sahara	closed	–	VL333SA	VL334SA	VL335SA	VL336SA	VL337SA
	black	closed	–	VL333BL	VL334BL	VL335BL	VL336BL	VL337BL
ATU Maternity Pantyhose	sahara	closed	–	VL343SA	VL344SA	VL345SA	VL346SA	VL347SA
	black	closed	–	VL343BL	VL344BL	VL345BL	VL346BL	VL347BL

VENOSAN® LEGLINE 30

ELEGANT LUXURY

Click for
VENOSAN
LEGLINE 30
product page




Return to
contents page

off-the-shelf range

Sizing Guide

Circumference	X-Small (cm)	Small (cm)	Medium (cm)	Large (cm)	X-Large (cm)
cB - Ankle	20 – 21	22 – 23	24 – 25	26 – 27	28 – 29
cC - Calf	31 – 33	34 – 36	37 – 39	40 – 42	43 – 45
cD - Below knee	29 – 31	32 – 34	35 – 37	38 – 40	41 – 43
cF - Mid thigh	42 – 44	45 – 47	48 – 52	51 – 53	54 – 60



Ordering Information

VENOSAN® LEGLINE 30

ELEGANT LUXURY

30mmHg

Style	Colour	Toe	Silicone Band	REF				
				XS	S	M	L	XL
AD Below knee	nude	closed	–	VL403NU	VL404NU	VL405NU	VL406NU	
	sahara	closed	–	VL403SA	VL404SA	VL405SA	VL406SA	
	black	closed	–	VL403BL	VL404BL	VL405BL	VL406BL	
AGH Thigh length with silicone	nude	closed	lace top	VL423NU	VL424NU	VL425NU	VL426NU	
	sahara	closed	lace top	VL423SA	VL424SA	VL425SA	VL426SA	
	black	closed	lace top	VL423BL	VL424BL	VL425BL	VL426BL	
AT Pantyhose	nude	closed	–	VL433NU	VL434NU	VL435NU	VL436NU	VL437NU
	sahara	closed	–	VL433SA	VL434SA	VL435SA	VL436SA	VL437SA
	black	closed	–	VL433BL	VL434BL	VL435BL	VL436BL	VL437BL



Return to contents page

VENOSAN® MICROFIBERLINE

STYLISH COMFORT, FASHIONABLY SOFT

off-the-shelf range

Sizing Guide

VENOSAN® MicroFiberLine Gentlemen

Circumference	Small (cm)	Medium (cm)	Large (cm)	X-Large (cm)
cB - Ankle	19 – 21.5	21.5 – 24	24 – 26.5	26.5 – 32
cC - Calf	28 – 38	33 – 43.5	35.5 – 46	38.5 – 51.5

Click for VENOSAN MICROFIBERLINE product page

VENOSAN® MicroFiberLine Lady

Circumference	Small (cm)	Medium (cm)	Large (cm)	X-Large (cm)
cB - Ankle	16.5 – 20.5	20.5 – 24.5	24.5 – 28	28 – 32
cC - Calf	26.5 – 37	29 – 39.5	32 – 43	34.5 – 47

Shoe size chart

VENOSAN® MicroFiberLine Gentlemen

	Small	Medium	Large	X-Large
US size	7 & under	7½ – 10	10½ – 12	12½ +

VENOSAN® MicroFiberLine Lady

	Small	Medium	Large	X-Large
US size	4 – 5	5½ – 7½	8 – 10½	11 – 13

Ordering Information

VENOSAN® MicroFiberLine Gentlemen

Style	Compression	Colour	Toe	REF			
				S	M	L	XL
AD Below knee	15 – 20mmHg	black	closed	MC79002	MC79003	MC79004	MC79005
	20 – 30mmHg	black	closed	MC89002	MC89003	MC89004	MC89005

VENOSAN® MicroFiberLine Lady

Style	Compression	Colour	Toe	REF			
				S	M	L	XL
AD Below knee	15 – 20mmHg	black	closed	MC79221	MC79222	MC79223	MC79224
		beige	closed	MC79205	MC79206	MC79207	MC79208

VENOSAN® SILVERLINE EVERYDAY COMFORT

Click for
VENOSAN
SILVERLINE
product page



Return to
contents page

off-the-shelf range

Sizing Guide

VENOSAN® Silverline Gentlemen

Circumference	Small (cm)	Medium (cm)	Large (cm)	X-Large (cm)	XX-Large (cm)
cB - Ankle	19 – 21.5	21.5 – 24	24 – 26.5	26.5 – 32	32 – 36
cC - Calf	28 – 38	33 – 43.5	35.5 – 46	38.5 – 51.5	41 – 53

VENOSAN® Silverline Lady

Circumference	Small (cm)	Medium (cm)	Large (cm)	X-Large (cm)	XX-Large (cm)
cB - Ankle	16.5 – 20.5	20.5 – 24.5	24.5 – 28	28 – 32	32 – 36
cC - Calf	26.5 – 37	29 – 39.5	32 – 43	34.5 – 47	41 – 53

Shoe size chart

VENOSAN® Silverline Gentlemen

	Small	Medium	Large	X-Large	XX-Large
US size	7 & under	7½ – 10	10½ – 12	12½ – 14	13 – 15

VENOSAN® Silverline Lady

	Small	Medium	Large	X-Large	XX-Large
US size	4 – 5	5½ – 7½	8 – 10½	11 – 13	14 – 15

Ordering Information

VENOSAN® Silverline Gentlemen

Style	Compression	Colour	Toe	REF				
				S	M	L	XL	XXL
AD Below knee	20-30mmHg	black	closed	VS49405	VS49406	VS49407	VS49408	VS49409

VENOSAN® Silverline Lady

Style	Compression	Colour	Toe	REF				
				S	M	L	XL	XXL
AD Below knee	20-30mmHg	black	closed	VS49505	VS49506	VS49507	VS49508	VS49509
		beige	closed	VS49511	VS49512	VS49513	VS49514	VS49515

VENOSAN® ATHLETIC PERFORMANCE SUPPORT

Click for
VENOSAN
ATHLETIC
product page



Return to
contents page

off-the-shelf range

Sizing Guide

Circumference	Small (cm)	Medium (cm)	Large (cm)	X-Large (cm)
cB - Ankle	18 – 21	22 – 25	25 – 29	29 – 33
cC - Calf	29 – 38	30 – 42	32 – 46	34 – 50

Shoe size chart

	Small	Medium	Large	X-Large
Mens' US	7 & under	7½ – 10	10½ – 12	12½+
Ladies' US	6 – 8	8½ – 11	11½ – 13	13+

Ordering Information

Style	Compression	Colour	Toe	REF			
				S	M	L	XL
AD Below knee	20 – 30mmHg	white	closed	A59202W	A59203W	A59204W	A59205W
		black	closed	A59202B	A59203B	A59204B	A59205B
		pink	closed	A59202F	A59203F	A59204F	A59205F
		purple	closed	A59202P	A59203P	A59204P	A59205P
		blue	closed	A59202S	A59203S	A59204S	A59205S

VENOSAN® SPORTSLINE PERFORMANCE SUPPORT

Click for
VENOSAN
SPORTSLINE
product page



Return to
contents page

off-the-shelf range

Sizing Guide

Circumference	Small (cm)	Medium (cm)	Large (cm)	X-Large (cm)
cB - Ankle	18 – 21	22 – 25	25 – 29	29 – 33
cC - Calf	29 – 38	30 – 42	32 – 46	34 – 50

Shoe size chart

	Small	Medium	Large	X-Large
Mens' US	7 & under	7½ – 10	10½ – 12	12½+
Ladies' US	6 – 8	8½ – 11	11½ – 13	13+

Ordering Information

Style	Compression	Colour	Toe	REF			
				S	M	L	XL
AD Below knee	15-20mmHg	white	closed	SP69202	SP69203	SP69204	SP69205
		black	closed	SP89202	SP89203	SP89204	SP89205
	20-30mmHg	white	closed	SP79202	SP79203	SP79204	SP79205
		black	closed	SP99202	SP99203	SP99204	SP99205

VENOSAN® AES

THROMBOSIS AND EMBOLISM PREVENTION

Click for
VENOSAN AES
product page



Return to
contents page

Ordering Information

18 mm Hg

Style	Ankle circumference cB	Calf circumference cC	Thigh circumference cG	Length	Size	Color coding Size (toe)	Color coding Length (top)	REF	
AD Below knee	18 – 20 cm	< 32 cm		< 41 cm	S normal	yellow	natural	SG57000	
				> 41 cm	S long	yellow	blue	SG57001	
	20 – 22 cm	32 – 37 cm		< 43 cm	M normal	natural	natural	SG57002	
				> 43 cm	M long	natural	blue	SG57003	
	23 – 25 cm	37 – 42 cm		< 45 cm	L normal	blue	natural	SG57004	
				> 45 cm	L long	blue	blue	SG57005	
	25 – 27 cm	42 – 47 cm		< 45 cm	XL normal	green	natural	SG57006	
				> 45 cm	XL long	green	blue	SG57007	
	28 – 30 cm	47 – 52 cm		< 45 cm	XXL normal	purple	natural	SG57008	
				> 45 cm	XXL long	purple	blue	SG57009	
	AGH Thigh length with silicone	18 – 20 cm	< 32 cm	45 – 64 cm	< 74 cm	S short	yellow	yellow	SG57020
					74 – 79 cm	S normal	yellow	natural	SG57021
> 79 cm					S long	yellow	blue	SG57022	
20 – 22 cm		32 – 37 cm	53 – 71 cm	< 74 cm	M short	natural	yellow	SG57023	
				74 – 79 cm	M normal	natural	natural	SG57024	
				> 79 cm	M long	natural	blue	SG57025	
23 – 25 cm		37 – 42 cm	61 – 80 cm	< 74 cm	L short	blue	yellow	SG57026	
				74 – 79 cm	L normal	blue	natural	SG57027	
				> 79 cm	L long	blue	blue	SG57028	
25 – 27 cm		42 – 47 cm	70 – 84 cm	< 79 cm	XL normal	green	natural	SG57029	
				> 79 cm	XL long	green	blue	SG57030	
AGG Thigh length with belt		18 – 20 cm	< 32 cm	45 – 64 cm	< 74 cm	S short	yellow	yellow	SG57050
	74 – 79 cm				S normal	yellow	natural	SG57051	
	> 79 cm				S long	yellow	blue	SG57052	
	20 – 22 cm	32 – 37 cm	53 – 71 cm	< 74 cm	M short	natural	yellow	SG57053	
				74 – 79 cm	M normal	natural	natural	SG57054	
				> 79 cm	M Long	natural	blue	SG57055	
	23 – 25 cm	37 – 42 cm	61 – 80 cm	< 74 cm	L short	blue	yellow	SG57056	
				74 – 79 cm	L normal	blue	natural	SG57057	
				> 79 cm	L long	blue	blue	SG57058	
	25 – 27 cm	42 – 47 cm	70 – 84 cm	< 79 cm	XL normal	green	natural	SG57059	
				> 79 cm	XL long	green	blue	SG57060	

ReadyWrap™, TributeWrap™, Caresia™, Exo Measuring Instructions



Return to
contents page


Lower Extremity

Circumference points:

- cG** 2 cm below gluteal fold
- cE¹** 5 cm above the knee
- cD** 2-3 cm below popliteal fossa
- cC** widest part of calf
- cB** 2 cm above malleolus
- cA¹** around the arch/dorsum of the foot
- cA** around base of toes

Lengths:

- IA-Y** heel to base of toes
- ID-B** the length between point cB and point cD
- IG-E¹** the length between point cE¹ (5 cm above the knee) and point cG




Upper Extremity

Circumference points:

- cA** widest width of the palm, just below the metacarpal phalangeal joints
- cC** narrowest point of the wrist
- cE** round the elbow joint and into the crease of the elbow
- cG** top of the arm, just below the axillary crease

Lengths:

- IA-C** the length between point cA and point cC
- IC-G** the length between point cC and point cG



ReadyWrap™

Upper Extremity Garments

Click for
ReadyWrap
product page


Return to
contents page

Ordering Information

ReadyWrap Arm

	Circumference (cm)			Length (cm) (IC-G)	REF	
	Wrist (cC)	Elbow (cE)	Axila (cG)		Beige	Black
Small	13.5 – 17	19 – 26	22 – 29	Short (up to 43cm)	104514	104515
				Average (up to 46cm)	104516	104517
				Tall (up to 49.5cm)	104518	104519
Medium	15.5 – 20	24 – 33	26 – 36	Short (up to 43cm)	104508	104509
				Average (up to 46cm)	104510	104511
				Tall (up to 49.5cm)	104512	104513
Large	17.5 – 23	29 – 40	31 – 43	Short (up to 43cm)	104502	104503
				Average (up to 46cm)	104504	104505
				Tall (up to 49.5cm)	104506	104507
X-Large	17.5 – 23	30 – 43	36 – 50	Short (up to 43cm)	104520	104521
				Average (up to 46cm)	104522	104523
				Tall (up to 49.5cm)	104524	104525



ReadyWrap Gauntlet

	Circumference (cm)		Length (cm) (IA – C)		REF	
	Palm (cA)	Wrist (cC)			Beige	Black
Small	18 – 20	13.5 – 17.0	up to 12.5	Right	104498	104499
				Left	104490	104491
Medium	20 – 22	15.5 – 20.0	up to 13.5	Right	104496	104497
				Left	104488	104489
Large	22 – 24	17.5 – 23.0	up to 14.5	Right	104494	104495
				Left	104486	104487
X-Large	24 – 26	19.5 – 26.0	up to 15.5	Right	104500	104501
				Left	104492	104493



ReadyWrap™ Lower Extremity Garments

Click for
ReadyWrap
product page



Return to
contents page

Ordering Information

ReadyWrap Toes

	Circumference (cm)			REF	
	Base of toes (cA)			Beige	Black
	Small	18.0 – 20.4		Right	104416
		Left	104406	104407	
Medium	20.5 – 23.5	Right	104414	104415	
		Left	104404	104405	
Large	23.6 – 26.0	Right	104412	104413	
		Left	104402	104403	
X-Large	26.1 – 28.5	Right	104418	104419	
		Left	104408	104409	
XX-Large	28.6 – 31.0	Right	104420	104421	
		Left	104410	104411	



ReadyWrap Foot CT

	Circumference (cm)			Foot length (cm) (A-Y)		REF	
	Base of toes (cA)	Foot (cA¹)	Ankle (cB)			Beige	Black
	Small	18.0-20.4	20-25			18-23	16-19cm (Regular)
					Left	141726	141728
				19-22cm (Long)	Right	141747	141749
					Left	141746	141748
Medium	20.5-23.5	24-29	20-28	16-19cm (Regular)	Right	141731	141733
					Left	141730	141732
				19-22cm (Long)	Right	141751	141753
					Left	141750	141752
Large	23.6-26.0	28-33	25-33	16-19cm (Regular)	Right	141735	141737
					Left	141734	141736
				19-22cm (Long)	Right	141755	141757
					Left	141754	141756
X-Large	26.1-28.5	32-37	30-38	16-19cm (Regular)	Right	141739	141741
					Left	141738	141740
				19-22cm (Long)	Right	141759	141761
					Left	141758	141760
XX-Large	28.6-31.0	36-41	35-43	16-19cm (Regular)	Right	141743	141745
					Left	141742	141744
				19-22cm (Long)	Right	141763	141765
					Left	141762	141764



ReadyWrap™

Lower Extremity Garments

Click for
ReadyWrap
product page



Return to
contents page

Ordering Information

ReadyWrap Foot SL

	Circumference (cm)		Foot length (cm) (IA-Y)		REF	
	Foot (cA ¹)	Ankle (cB)			Beige	Black
Small	21 – 24	18 – 23	16 – 19cm (Regular)	Right	144425	144427
				Left	144424	144426
			19 – 22cm (Long)	Right	144429	144431
				Left	144428	144430
Medium	25 – 28	20 – 28	16 – 19cm (Regular)	Right	144433	144435
				Left	144432	144434
			19 – 22cm (Long)	Right	144437	144439
				Left	144436	144438
Large	29 – 32	25 – 33	16 – 19cm (Regular)	Right	144441	144443
				Left	144440	144442
			19 – 22cm (Long)	Right	144445	144447
				Left	144444	144446
X-Large	33 – 36	30 – 38	16 – 19cm (Regular)	Right	144449	144451
				Left	144448	144450
			19 – 22cm (Long)	Right	144453	144455
				Left	144452	144454



ReadyWrap Calf

	Circumference (cm)		Length (cm) (IB-D)	REF	
	Ankle (cB)	Calf (cC)		Beige	Black
Small	18 – 23	28 – 38	up to 30	104450	104451
			up to 35	104452	104453
Medium	20 – 28	33 – 43	up to 30	104446	104447
			up to 35	104448	104449
Large	25 – 33	38 – 48	up to 30	104442	104443
			up to 35	104444	104445
X-Large	30 – 38	43 – 58	up to 30	104454	104455
			up to 35	104456	104457
XX-Large	35 – 43	53 – 68	up to 30	104784	104785
			up to 35	104786	104787



ReadyWrap Knee

	Circumference (cm)		REF	
	Below Knee (cD)	Above Knee (cE ¹)	Beige	Black
Small	28 – 38	38 – 48	104464	104465
Medium	33 – 43	45 – 58	104462	104463
Large	38 – 48	52 – 68	104460	104461
X-Large	43 – 58	59 – 78	104466	104467
XX-Large	53 – 68	66 – 88	104788	104789



ReadyWrap™ Lower Extremity Garments

Click for
ReadyWrap
product page



Return to
contents page

Ordering Information

ReadyWrap Thigh

	Circumference (cm)		Length (cm) (IE' - G)	REF	
	Above Knee (cE')	Groin (cG)		Beige	Black
Small	38-48	48-63	up to 25	104478	104479
			up to 30	104480	104481
Medium	45-58	58-73	up to 25	104474	104475
			up to 30	104476	104477
Large	52-68	68-83	up to 25	104470	104471
			up to 30	104472	104473
X-Large	59-78	78-93	up to 25	104482	104483
			up to 30	104484	104485
XX-Large	66-88	88-103	up to 25	104790	104791
			up to 30	104792	104793



ReadyWrap Liners

For ReadyWrap below knee garments

	REF
	Black
Small – Large	149032
X-Large –XX-Large	149033

For ReadyWrap Thigh garments

	REF
	Black
Small – Large	149034
X-Large –XX-Large	149035



ReadyWrap Extender Strap

	REF	
	Beige	Black
One Size	304398	304399



TributeWrap™

Upper & Lower Extremity Garments

Ordering Information

TributeWrap Glove

Size	Circumference (cm)			REF
	Palm (cB)	Wrist (cC)		
S	16 – 20	16.5 – 21	Left	146783
			Right	146786
M	18 – 22	17.5 – 22	Left	146784
			Right	146787
L	20 – 24.5	18.5 – 23	Left	146785
			Right	146788



TributeWrap Arm

Size	Circumference (cm)			Length (cm)	REF	
	Wrist (cC)	Elbow (cE)	Axilla (cG)			
S	16 – 21	22 – 28	25 – 31	43 – 47	Left	146777
					Right	146780
M	17 – 22	26 – 32	30 – 37	43 – 47	Left	146778
					Right	146781
L	18 – 23	30 – 36	35 – 44	43 – 47	Left	146779
					Right	146782



TributeWrap Below Knee

Size	Circumference (cm)		Length (cm)		REF	
	Ankle (cB)	Calf (cC)	Foot (IA)	Calf (ID)		
S	24 – 30	30 – 36	19 – 22	40 – 44	Left	146789
					Right	146792
M	27 – 33	33 – 41.5	19 – 22	40 – 44	Left	146790
					Right	146793
L	30 – 36	36 – 45	19 – 22	40 – 44	Left	146791
					Right	146794



Caresia™

Upper Extremity Garments

Click for
Caresia
product page



Return to
contents page

Ordering Information

Palm to Axilla

Description	Arm Length: 59 – 65cm (IC-G)			REF
	Palm Circumference (cA)	Wrist Circumference (cC)	Biceps Circumference (cF)	
Small (Left)	15 – 19cm	13.5 – 17cm	20 – 30cm	104384
Small (Right)	15 – 19cm	13.5 – 17cm	20 – 30cm	104387
Medium (Left)	19 – 23cm	17 – 22cm	30 – 40cm	104383
Medium (Right)	19 – 23cm	17 – 22cm	30 – 40cm	104386
Large (Left)	23 – 28cm	22 – 28cm	40 – 50cm	104382
Large (Right)	23 – 28cm	22 – 28cm	40 – 50cm	104385



MCP to Axilla


Wrist to Axilla

Description	Arm Length: 47–53cm (IC-G)		REF
	Wrist Circumference (cC)	Biceps Circumference (cF)	
Small (Left)	13.5 – 17cm	20 – 30cm	104390
Small (Right)	13.5 – 17cm	20 – 30cm	104393
Medium (Left)	17 – 22cm	30 – 40cm	104389
Medium (Right)	17 – 22cm	30 – 40cm	104392
Large (Left)	22 – 28cm	40 – 50cm	104388
Large (Right)	22 – 28cm	40 – 50cm	104391



Glove

Description	Length: up to 26cm (IA-C)		REF
	Palm Circumference (cA)	Wrist Circumference (cC)	
Small	15 – 19cm	13.5 – 17cm	104378
Medium	19 – 23cm	17 – 22cm	104377
Large	23 – 28cm	22 – 28cm	104376



Gauntlet

Description	Length: up to 18cm (IA-C)		REF
	Palm Circumference (cA)	Wrist Circumference (cC)	
Small	15 – 19cm	13.5 – 17cm	104381
Medium	19 – 23cm	17 – 22cm	104380
Large	23 – 28cm	22 – 28cm	104379



Caresia™

Lower Extremity Garments

Click for
Caresia
product page



Return to
contents page

Ordering Information

Thigh

Short	Thigh Length: 28 – 36cm (IG-E1)		REF
	Groin Circumference (cG)		
Small	50 – 70cm		104370
Medium	60 – 80cm		104369
Large	70 – 90cm		104368
Standard	Thigh Length: 33 – 41cm (IG-E1)		REF
	Groin Circumference (cG)		
Small	50 – 70cm		104366
Medium	60 – 80cm		104365
Large	70 – 90cm		104364
Tall	Thigh Length: 39 – 47cm (IG-E1)		REF
	Groin Circumference (cG)		
Small	50 – 70cm		104374
Medium	60 – 80cm		104373
Large	70 – 90cm		104372



Below Knee

Short	Calf Length: 34 – 40cm (ID)		REF
	Ankle Circumference (cB)	Calf Circumference (cC)	
Small	22 – 28cm	30 – 40cm	104360
Medium	28 – 34cm	40 – 50cm	104359
Large	34 – 40cm	50 – 60cm	104358
Standard	Calf Length: 40 – 46cm (ID)		REF
	Ankle Circumference (cB)	Calf Circumference (cC)	
Small	22 – 28cm	30 – 40cm	104357
Medium	28 – 34cm	40 – 50cm	104356
Large	34 – 40cm	50 – 60cm	104355
Tall	Calf Length: 46 – 52cm (ID)		REF
	Ankle Circumference (cB)	Calf Circumference (cC)	
Small	22 – 28cm	30 – 40cm	104363
Medium	28 – 34cm	40 – 50cm	104362
Large	34 – 40cm	50 – 60cm	104361



Foot

Description	Foot Length (IZ)	Ankle Circumference (cB)	REF
Small	20 – 25cm	22 – 28cm	104354
Medium	25 – 30cm	28 – 34cm	104353
Large	30 – 35cm	34 – 40cm	104352



Ordering Information

Image Reference	Product Name	Width (cm)	Height (cm)	Length (cm)	Diameter (cm)	REF
A	Breast Pad (A-B cup)	-	-	-	-	104541
A	Breast Pad (C-D cup)	-	-	-	-	104540
A	Breast Pad (D-DD cup)	-	-	-	-	104539
B	Carpal Hand Pad	4.5	7	-	-	104537
C	Circular Pad Small	-	-	-	18	104536
C	Circular Pad Large	-	-	-	22	104535
D	Crescent Pad Small	-	6.5	51	9	104547
D	Crescent Pad Medium	-	6.5	61	11	104546
D	Crescent Pad Large	-	6.5	71	12	104545
	Curved Pad	5	13	-	-	104527
E	Dorsal Pad Small	10	14	-	-	104529
E	Dorsal Pad Large	10	19	-	-	104528
F	Female Genital Pad	7	5	-	-	104560
	Full Bra Pad (A-B cup)	-	-	-	-	104556
	Full Bra Pad (C-D cup)	-	-	-	-	104555
	Full Bra Pad (C-D cup)	-	-	-	-	104554
G	Head-Neck Pad	16.5	7.5	-	-	104544
H	Kidney Pad	14	24	-	-	104538
I	Knee-Elbow Shaper Pad	Girth range 35 - 65			-	104557
J	Labia Pad	33	53	-	-	104551
K	Lateral Bra Pad	10	16	-	-	104553
L	Oval Pad	8	19	-	-	104532
M	Palmar Pad Small	7	8	-	-	104531
M	Palmar Pad Large	8	9	-	-	104530
N	Rectangular-Tapered Pad Small	11	36	-	-	104534
N	Rectangular-Tapered Pad Large	18	42	-	-	104533
O	Ring Pad (with VELCRO)	60	6.5	-	9	104558
O	Ring Pad (without VELCRO)	60	6.5	-	9	104559
	Scrotum Pad	34	53	-	-	104552
	Thorax-Thigh Pad	52	22	-	-	104543
	Tic-Tac Pad Small	20	-	20	-	104563
	Tic-Tac Pad Medium	30	-	30	-	104562
	Tic-Tac Pad Large	40	-	40	-	104561
P	Toe Pad Small	MTP circumference <22		-	-	104550
P	Toe Pad Medium	MTP circumference <25		-	-	104549
P	Toe Pad Large	MTP circumference <32		-	-	104584
Q	Versi-Face Pad	29	16	-	-	104542



Ordering Information

ExoSoft Glove

sold individually

Compression 20-30mmHg

Size	Finger Length							
	Circumference (cm)		Quarter		Half		Full	
	Palm (cB)	Wrist (cC)	Beige	Black	Beige	Black	Beige	Black
	REF							
Small	17 – 19	15 – 17	104116	104119	104117	104120	104118	104121
Medium	19 – 21	17 – 19	104122	104125	104123	104126	104124	104127
Large	22 – 24	20 – 22	104128	104131	104129	104132	104130	104133
X-Large	25 – 27	23 – 25	104134	104138	104135	104139	104136	104140



ExoSoft Gauntlet

sold individually

Compression 20-30mmHg

	Circumference (cm)		REF	
	Palm (cB)	Wrist (cC)	Beige	Black
Small	17 – 19	15 – 17	104141	104142
Medium	19 – 21	17 – 19	104143	104144
Large	22 – 24	20 – 22	104145	104146
X-Large	25 – 27	23 – 25	104147	104148



ExoSoft Armsleeve

sold individually

Compression 20-30mmHg

Size	Length (cm)	Top of Armsleeve											
		Circumference (cm)			Knit		Silicone						
		(IG)	Wrist (cC)	Elbow (cE)	Axilla (cG)	Beige	Black	Beige	Black				
		REF											
Small	Short (up to 43cm)	15 – 17	24 – 27	27 – 30	104149	104151	104150	104152					
	Tall (more than 43cm)				104165	104167	104166	104168					
Medium	Short (up to 43cm)	17 – 19	27 – 31	31 – 35	104153	104155	104154	104156					
	Tall (more than 43cm)				104169	104171	104170	104172					
Large	Short (up to 43cm)	20 – 22	31 – 35	35 – 40	104157	104159	104158	104160					
	Tall (more than 43cm)				104173	104175	104174	104176					
X-Large	Short (up to 43cm)	23 – 25	36 – 41	40 – 45	104161	104163	104162	104164					
	Tall (more than 43cm)				104177	104178	104179	104180					



ExoStrong™ Upper Extremity Garments

Click for
ExoStrong
product page



Return to
contents page

Ordering Information

ExoStrong Gauntlet

sold individually

Compression 20-30mmHg

	Circumference (cm)		REF	
	Palm (cB)	Wrist (cC)	Beige	Black
	Small	17 – 19	15 – 17	104296
Medium	19 – 21	17 – 19	104298	104299
Large	22 – 24	20 – 22	104300	104301
X-Large	25 – 27	23 – 25	104302	104303



ExoStrong Glove

sold individually

Compression 20-30mmHg

Size	Circumference (cm)		Finger Length					
			Quarter		Half		Full	
	Palm (cB)	Wrist (cC)	Beige	Black	Beige	Black	Beige	Black
			REF					
Small	17 – 19	15 – 17	104272	104275	104273	104276	104274	104277
Medium	19 – 21	17 – 19	104278	104281	104279	104282	104280	104283
Large	22 – 24	20 – 22	104284	104287	104285	104288	104286	104289
X-Large	25 – 27	23 – 25	104290	104293	104291	104294	104292	104295



ExoStrong Armsleeve

sold individually

Compression 20-30mmHg

Size	Length (cm)	Top of Armsleeve						
		Circumference (cm)			Knit		Silicone	
	(IG)	Wrist (cC)	Elbow (cE)	Axilla (cG)	Beige	Black	Beige	Black
	REF							
Small	Short (up to 43cm)	15 – 17	24 – 27	27 – 30	104304	104305	104306	104307
	Average (up to 46cm)				104320	104322	104321	104323
	Tall (up to 49cm)				104336	104338	104337	104339
Medium	Short (up to 43cm)	17 – 19	27 – 31	31 – 35	104308	104310	104309	104311
	Average (up to 46cm)				104324	104326	104325	104327
	Tall (up to 49cm)				104340	104342	104341	104343
Large	Short (up to 43cm)	20 – 22	31 – 35	35 – 40	104312	104314	104313	104315
	Average (up to 46cm)				104328	104330	104328	104331
	Tall (up to 49cm)				104344	104346	104345	104347
X-Large	Short (up to 43cm)	23 – 25	36 – 41	40 – 45	104316	104318	104317	104319
	Average (up to 46cm)				104332	104334	104333	104335
	Tall (up to 49cm)				104348	104350	104349	104351



ExoStrong™ Lower Extremity Garments

Click for
ExoStrong
product page



Return to
contents page

Ordering Information

ExoStrong Below Knee

sold individually

Compression 20-30mmHg Top Style: Open

Size	Length (cm)		Circumference (cm)		Top of Stocking				
					Knit		Silicone		
	(IA)	(ID)	(cB)	(cC)	Beige	Black	Beige	Black	
									REF
Small	Short	up to 17	up to 38	18 – 20	27 – 36	104181	104182	104183	104184
	Average	up to 19	up to 41			104201	104203	104202	104204
	Tall	up to 22	up to 44			104221	104223	104222	104224
Medium	Short	up to 17	up to 38	20 – 24	29 – 39	104185	104187	104186	104188
	Average	up to 19	up to 41			104205	104207	104206	104208
	Tall	up to 22	up to 44			104225	104227	104226	104228
Large	Short	up to 17	up to 38	24 – 28	32 – 43	104189	104191	104190	104192
	Average	up to 19	up to 41			104209	104211	104210	104212
	Tall	up to 22	up to 44			104229	104231	104230	104232
X-Large	Short	up to 17	up to 38	28 – 31	33 – 45	104193	104195	104194	104196
	Average	up to 19	up to 41			104213	104215	104214	104216
	Tall	up to 22	up to 44			104233	104235	104234	104236
XX-Large	Short	up to 17	up to 38	32 – 34	37 – 50	104197	104199	104198	104200
	Average	up to 19	up to 41			104217	104219	104218	104220
	Tall	up to 22	up to 44			104237	104239	104238	104240



ExoStrong Thigh High

sold individually

Compression 20-30mmHg Top Style: Open

Size	Length (cm)		Circumference (cm)		Top of Stocking			
					Silicone			
	(IA)	(ID)	(cB)	(cC)	(cG)	Beige	Black	
								REF
Small	Short	up to 17	up to 71	18 – 20	27 – 36	41 – 55	104241	104242
	Average	up to 19	up to 74				104251	104252
	Tall	up to 22	up to 77				104261	104262
Medium	Short	up to 17	up to 71	20 – 24	29 – 39	45 – 61	104243	104244
	Average	up to 19	up to 74				104253	104254
	Tall	up to 22	up to 77				104263	104264
Large	Short	up to 17	up to 71	24 – 28	32 – 43	50 – 66	104245	104246
	Average	up to 19	up to 74				104255	104256
	Tall	up to 22	up to 77				104265	104266
X-Large	Short	up to 17	up to 71	28 – 31	33 – 45	56 – 71	104247	104248
	Average	up to 19	up to 74				104257	104258
	Tall	up to 22	up to 77				104267	104268
XX-Large	Short	up to 17	up to 71	32 – 34	37 – 50	66 – 81	104249	104250
	Average	up to 19	up to 74				104259	104260
	Tall	up to 22	up to 77				104269	104270



Measuring and ordering information



Lohmann & Rauscher

Rosidal® and ReadyWrap™ are trademarks of Lohmann & Rauscher GmbH & Co KG. © 2023

VENOSAN® is a registered trademark of Swisslastic AG St. Gallen

Lohmann & Rauscher Pty Ltd · 321 Ferntree Gully Road · Mount Waverley, Victoria 3149, Australia · **Customer Service:** 1300 572 869

Lohmann & Rauscher Limited · 21 Queen Street · Auckland 1010 · New Zealand · **Customer Service:** 0800 572 869