


- Single use - reduces the risk of cross infection.
- High gloss, low-friction outer surface improves tube glide over bougie.
- Low point-of-pressure tip reduces risk of accidental trauma.
- Highly visible colour for improved positioning in poor light.
- Extra tip graduations allow for improved laryngeal inlet entry and exit guidance.
- Adult size available in 600mm and 800mm length, each with or without angled (coudé) tip, allowing greater choice for clinicians. All available in 10ch and 15ch, suitable for use with all tubes from 4mm(ID) to 10mm(ID).
- Straight paediatric bougie, 5ch, 470mm length, also available.
- Individually packed sterile, 10 units within highly visible outer tube for ease of product identification in emergency.
- Robust outer tube prevents damage in transit and storage.

- Barium tip for improved X-Ray guidance if required (adult sizes only).



Code	Description	Size
S05470	Bougie, straight	5ch x 470mm
S10600	Bougie, straight	10ch x 600mm
S15600	Bougie, straight	15ch x 600mm
S10800	Bougie, straight	10ch x 800mm
S15800	Bougie, straight	15ch x 800mm
A10600	Bougie, angled (coudé)	10ch x 600mm
A15600	Bougie, angled (coudé)	15ch x 600mm
A10800	Bougie, angled (coudé)	10ch x 800mm
A15800	Bougie, angled (coudé)	15ch x 800mm

